

Shamir Product Family

To view a specific product material availability, select the product name from the list below.

To return to this page, click on the Shamir logo located in the bottom right-hand corner of the page you are on.

Everyday> Progressive

Shamir Autograph III®

Shamir Autograph II®

Shamir InTouch™

Shamir Spectrum™

Shamir Element™

Shamir FirstPAL™

Shamir Creation®

Shamir Piccolo®

Shamir Genesis®

Everyday> Bifocal

Shamir Duo™

Everyday> Single Vision

Shamir Autograph III® - SV

Shamir Autograph II - Single Vision™

Shamir Spectrum - Single Vision™

Shamir Single Vision™

Performance

Shamir Attitude III® - Fashion

Performance> Sport

Shamir Attitude III® - Sport

Shamir Golf™

Shamir Autograph II - Attitude™

Shamir Attitude®

Performance> Sport Single Vision

Shamir Attitude III® - SV

Shamir Golf - Single Vision™

Shamir Autograph II - Attitude Single Vision™

Occupation

Shamir Computer™ and WorkSpace™

Shamir Autograph II - Office™

Shamir Relax™

Shamir Office™



Shamir Autograph III®

PRESCRIPTOR ORDER CODE:

AUTOGRAPH III 11 AUTOGRAPH III 18
 AUTOGRAPH III 13 AUTOGRAPH III VARIABLE
 AUTOGRAPH III 15

PRODUCT FEATURES:

☒ FREEFORM PAL ☒ EYEPOINT TECHNOLOGY III® ☐ EYEPOINT TECHNOLOGY®
☐ FREEFORM SV ☒ AS-WORN QUADRO™ ☐ AS-WORN TECHNOLOGY™
☐ CONVENTIONAL PAL ☒ INTELLICORRIDOR™
☐ CONVENTIONAL SV ☒ NATURAL POSTURE™

NOTE:

- Recommended min fitting heights are 11mm, 13mm, 15mm and 18mm
- Variable 11mm and up
- Prism to 10.00 Diopters

CR-39	TREATMENT						COATING
	Clear	Transitions Signature7™	XTRActive	Vantage	Polarized	Drivewear	
POWER RANGE	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	Clear and Transitions
ADD POWER	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	
CYLINDER	to -6.00	to -6.00	to -6.00	to -6.00	to -6.00	to -6.00	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15		
INVISIBLE MARKINGS - 11mm	15 250 A ₁ ¹¹	TR 250 A ₁ ¹¹	X5 250 A ₁ ¹¹	TV 250 A ₁ ¹¹	PPZ 250 A ₁ ¹¹	DW 250 A ₁ ¹¹	
INVISIBLE MARKINGS - 13mm	15 250 A ₁ ¹³	TR 250 A ₁ ¹³	X5 250 A ₁ ¹³	TV 250 A ₁ ¹³	PPZ 250 A ₁ ¹³	DW 250 A ₁ ¹³	
INVISIBLE MARKINGS - 15mm	15 250 A ₁ ¹⁵	TR 250 A ₁ ¹⁵	X5 250 A ₁ ¹⁵	TV 250 A ₁ ¹⁵	PPZ 250 A ₁ ¹⁵	DW 250 A ₁ ¹⁵	Clear and Transitions
INVISIBLE MARKINGS - 18mm	15 250 A ₁ ¹⁸	TR 250 A ₁ ¹⁸	X5 250 A ₁ ¹⁸	TV 250 A ₁ ¹⁸	PPZ 250 A ₁ ¹⁸	DW 250 A ₁ ¹⁸	
INVISIBLE MARKINGS - Variable	15 250 A ₃	TR 250 A ₃	X5 250 A ₃	TV 250 A ₃	PPZ 250 A ₃	DW 250 A ₃	

SuperLite™ 1.60	TREATMENT				COATING
	Clear	Transitions Signature7™	XTRActive	Polarized	
POWER RANGE	-15.00 to +8.00	-15.00 to +8.00	-15.00 to +8.00	-15.00 to +8.00	Clear and Transitions
ADD POWER	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	
CYLINDER	to -7.00	to -7.00	to -7.00	to -7.00	
OTHER		Gray/Brown	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS - 11mm	16 250 A ₁ ¹¹	T6 250 A ₁ ¹¹	X6 250 A ₁ ¹¹	P26 250 A ₁ ¹¹	
INVISIBLE MARKINGS - 13mm	16 250 A ₁ ¹³	T6 250 A ₁ ¹³	X6 250 A ₁ ¹³	P26 250 A ₁ ¹³	
INVISIBLE MARKINGS - 15mm	16 250 A ₁ ¹⁵	T6 250 A ₁ ¹⁵	X6 250 A ₁ ¹⁵	P26 250 A ₁ ¹⁵	Clear and Transitions
INVISIBLE MARKINGS - 18mm	16 250 A ₁ ¹⁸	T6 250 A ₁ ¹⁸	X6 250 A ₁ ¹⁸	P26 250 A ₁ ¹⁸	
INVISIBLE MARKINGS - Variable	16 250 A ₃	T6 250 A ₃	X6 250 A ₃	P26 250 A ₃	

SuperLite™ 1.67	TREATMENT				COATING
	Clear	Transitions Signature7™	XTRActive	Polarized	
POWER RANGE	-16.75 to +9.00	-16.75 to +9.00	-16.75 to +9.00	-16.75 to +9.00	Clear and Transitions
ADD POWER	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	
CYLINDER	to -8.00	to -8.00	to -8.00	to -8.00	
OTHER		Gray/Brown	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS - 11mm	67 250 A ₁ ¹¹	T67 250 A ₁ ¹¹	X67 250 A ₁ ¹¹	P267 250 A ₁ ¹¹	
INVISIBLE MARKINGS - 13mm	67 250 A ₁ ¹³	T67 250 A ₁ ¹³	X67 250 A ₁ ¹³	P267 250 A ₁ ¹³	
INVISIBLE MARKINGS - 15mm	67 250 A ₁ ¹⁵	T67 250 A ₁ ¹⁵	X67 250 A ₁ ¹⁵	P267 250 A ₁ ¹⁵	Clear and Transitions
INVISIBLE MARKINGS - 18mm	67 250 A ₁ ¹⁸	T67 250 A ₁ ¹⁸	X67 250 A ₁ ¹⁸	P267 250 A ₁ ¹⁸	
INVISIBLE MARKINGS - Variable	67 250 A ₃	T67 250 A ₃	X67 250 A ₃	P267 250 A ₃	

SuperLite™ 1.74	TREATMENT		COATING
	Clear	Transitions Signature7™	
POWER RANGE	-18.75 to +15.50	-18.75 to +15.50	Clear and Transitions
ADD POWER	+0.75 to +4.00	+0.75 to +4.00	
CYLINDER	to -9.00	to -9.00	
OTHER		Gray/Brown	
INVISIBLE MARKINGS - 11mm	74 250 A ₁ ¹¹	T74 250 A ₁ ¹¹	
INVISIBLE MARKINGS - 13mm	74 250 A ₁ ¹³	T74 250 A ₁ ¹³	
INVISIBLE MARKINGS - 15mm	74 250 A ₁ ¹⁵	T74 250 A ₁ ¹⁵	Clear and Transitions
INVISIBLE MARKINGS - 18mm	74 250 A ₁ ¹⁸	T74 250 A ₁ ¹⁸	
INVISIBLE MARKINGS - Variable	74 250 A ₃	T74 250 A ₃	

DLC™/TRIVEX™	TREATMENT					COATING
	Clear	Transitions Signature7™	XTRActive	Vantage	Polarized	
POWER RANGE	-13.25 to +7.00	-13.25 to +7.00	-13.25 to +7.00	-13.25 to +7.00	-13.25 to +7.00	Clear and Transitions
ADD POWER	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	
CYLINDER	to -6.50	to -6.50	to -6.50	to -6.50	to -6.50	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS - 11mm	TX 250 A ₁ ¹¹	TTX 250 A ₁ ¹¹	XX 250 A ₁ ¹¹	VX 250 A ₁ ¹¹	PTX 250 A ₁ ¹¹	
INVISIBLE MARKINGS - 13mm	TX 250 A ₁ ¹³	TTX 250 A ₁ ¹³	XX 250 A ₁ ¹³	VX 250 A ₁ ¹³	PTX 250 A ₁ ¹³	
INVISIBLE MARKINGS - 15mm	TX 250 A ₁ ¹⁵	TTX 250 A ₁ ¹⁵	XX 250 A ₁ ¹⁵	VX 250 A ₁ ¹⁵	PTX 250 A ₁ ¹⁵	Clear and Transitions
INVISIBLE MARKINGS - 18mm	TX 250 A ₁ ¹⁸	TTX 250 A ₁ ¹⁸	XX 250 A ₁ ¹⁸	VX 250 A ₁ ¹⁸	PTX 250 A ₁ ¹⁸	
INVISIBLE MARKINGS - Variable	TX 250 A ₃	TTX 250 A ₃	XX 250 A ₃	VX 250 A ₃	PTX 250 A ₃	

1.56	TREATMENT		COATING
	Clear	Polarized	
POWER RANGE	-14.00 to +6.50	-14.00 to +6.50	Clear and Transitions
ADD POWER	+0.75 to +4.00	+0.75 to +4.00	
CYLINDER	to -6.75	to -6.75	
OTHER	BluTech Indoor	BluTech Outdoor	
INVISIBLE MARKINGS - 11mm	56 250 A ₁ ¹¹	56P 250 A ₁ ¹¹	
INVISIBLE MARKINGS - 13mm	56 250 A ₁ ¹³	56P 250 A ₁ ¹³	
INVISIBLE MARKINGS - 15mm	56 250 A ₁ ¹⁵	56P 250 A ₁ ¹⁵	Clear and Transitions
INVISIBLE MARKINGS - 18mm	56 250 A ₁ ¹⁸	56P 250 A ₁ ¹⁸	
INVISIBLE MARKINGS - Variable	56 250 A ₃	56P 250 A ₃	

Crizal AR Coating Hoya AR Coating Zeiss AR Coating iCoat AR Coating

All Shamir products are available with your lab's in-house AR Coatings. Please check with your lab to determine which specific manufacturer's AR coating is available with the Shamir lens and material that you are ordering, as a manufacturer may offer more than one AR coating.

EFFECTIVE DATE: SEP. 15, 2014

- Denotes that material is available in a Shamir Freeform Blank™



SHAMIR

EVERYDAY> PROGRESSIVE

Shamir Autograph III®

PRESCRIPTOR ORDER CODE:

AUTOGRAPH III 11
AUTOGRAPH III 13
AUTOGRAPH III 15

AUTOGRAPH III 18
AUTOGRAPH III VARIABLE

PRODUCT FEATURES:

☒ FREEFORM PAL

☐ FREEFORM SV

☐ CONVENTIONAL PAL

☐ CONVENTIONAL SV

☒ EYEPOINT TECHNOLOGY III®

☒ AS-WORN QUADRO™

☒ INTELLICORRIDOR™

☒ NATURAL POSTURE™

☐ EYEPOINT TECHNOLOGY®

☐ AS-WORN TECHNOLOGY™

NOTE:

- Recommended min fitting heights are 11mm, 13mm, 15mm and 18mm
- Variable 11mm and up
- Prism to 10.00 Diopters

POLYPLUS™/ POLYCARBONATE	TREATMENT						COATING
	Clear	Transitions Signature7™	XTRActive	Vantage	Polarized	Drivewear	
POWER RANGE	-14.75 to +8.00	-14.75 to +8.00	-14.75 to +8.00	-14.75 to +8.00	-14.75 to +8.00	-14.75 to +8.00	<div>Clear and Transitions</div>
ADD POWER	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	
CYLINDER	to -7.00	to -7.00	to -7.00	to -7.00	to -7.00	to -7.00	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15		
INVISIBLE MARKINGS - 11mm	PC 250 A113	PT 250 A113	PX 250 A113	PVT 250 A113	PPZ 250 A113	DW 250 A113	
INVISIBLE MARKINGS - 13mm	PC 250 A133	PT 250 A133	PX 250 A133	PVT 250 A133	PPZ 250 A133	DW 250 A133	
INVISIBLE MARKINGS - 15mm	PC 250 A153	PT 250 A153	PX 250 A153	PVT 250 A153	PPZ 250 A153	DW 250 A153	
INVISIBLE MARKINGS - 18mm	PC 250 A183	PT 250 A183	PX 250 A183	PVT 250 A183	PPZ 250 A183	DW 250 A183	<div></div>
INVISIBLE MARKINGS - Variable	PC 250 A3	PT 250 A3	PX 250 A3	PVT 250 A3	PPZ 250 A3	DW 250 A3	

NXT	TREATMENT		COATING
	Photochromic	Polarized	
POWER RANGE	-1.50 to +7.00	-1.50 to +7.00	<div></div>
ADD POWER	+0.75 to +4.00	+0.75 to +4.00	
CYLINDER	to -6.50	to -6.50	
OTHER	Green/Copper	Gray/Brown	
INVISIBLE MARKINGS - 11mm	NXTP 250 A113	NXTZ 250 A113	
INVISIBLE MARKINGS - 13mm	NXTP 250 A133	NXTZ 250 A133	
INVISIBLE MARKINGS - 15mm	NXTP 250 A153	NXTZ 250 A153	
INVISIBLE MARKINGS - 18mm	NXTP 250 A183	NXTZ 250 A183	<div></div>
INVISIBLE MARKINGS - Variable	NXTP 250 A3	NXTZ 250 A3	

Crizal AR Coating Hoya AR Coating Zeiss AR Coating iCoat AR Coating

All Shamir products are available with your lab's in-house AR Coatings. Please check with your lab to determine which specific manufacturer's AR coating is available with the Shamir lens and material that you are ordering, as a manufacturer may offer more than one AR coating.

EFFECTIVE DATE: SEP. 15, 2014

- Denotes that material is available in a Shamir Freeform Blank™



SHAMIR
Autograph II®

PRESCRIPTOR ORDER CODE:

AUTOGRAPH II 11 AUTOGRAPH II 15
AUTOGRAPH II 13 AUTOGRAPH II 18

PRODUCT FEATURES:

- ☒ FREEFORM PAL
☐ FREEFORM SV
☐ CONVENTIONAL PAL
☐ CONVENTIONAL SV
- ☒ EYEPOINT TECHNOLOGY®
☒ AS-WORN TECHNOLOGY™

NOTE:

- Recommended min fitting heights are 11mm, 13mm, 15mm and 18mm
- Prism to 10.00 Diopters

CR-39	TREATMENT						COATING
	Clear	Transitions Signature7™	XTRActive	Vantage	Polarized	Drivewear	
POWER RANGE	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	
ADD POWER	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	
CYLINDER	to -6.00	to -6.00	to -6.00	to -6.00	to -6.00	to -6.00	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15		
INVISIBLE MARKINGS - 11mm	15 250 67	TR 250 67	X5 250 67	V5 250 67	PZ 250 67	DW 250 67	
INVISIBLE MARKINGS - 13mm	15 250 67	TR 250 67	X5 250 67	V5 250 67	PZ 250 67	DW 250 67	
INVISIBLE MARKINGS - 15mm	15 250 67	TR 250 67	X5 250 67	V5 250 67	PZ 250 67	DW 250 67	
INVISIBLE MARKINGS - 18mm	15 250 67	TR 250 67	X5 250 67	V5 250 67	PZ 250 67	DW 250 67	
INVISIBLE MARKINGS - Variable	15 250 67	TR 250 67	X5 250 67	V5 250 67	PZ 250 67	DW 250 67	

SuperLite™ 1.67	TREATMENT				COATING
	Clear	Transitions Signature7™	XTRActive	Polarized	
POWER RANGE	-16.75 to +8.25	-16.75 to +8.25	-16.75 to +8.25	-16.75 to +8.25	
ADD POWER	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	
CYLINDER	to -8.00	to -8.00	to -8.00	to -8.00	
OTHER		Gray/Brown	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS - 11mm	67 250 67	T67 250 67	X67 250 67	PZ67 250 67	
INVISIBLE MARKINGS - 13mm	67 250 67	T67 250 67	X67 250 67	PZ67 250 67	
INVISIBLE MARKINGS - 15mm	67 250 67	T67 250 67	X67 250 67	PZ67 250 67	
INVISIBLE MARKINGS - 18mm	67 250 67	T67 250 67	X67 250 67	PZ67 250 67	
INVISIBLE MARKINGS - Variable	67 250 67	T67 250 67	X67 250 67	PZ67 250 67	

SuperLite™ 1.60	TREATMENT				COATING
	Clear	Transitions Signature7™	XTRActive	Polarized	
POWER RANGE	-15.00 to +6.00	-15.00 to +6.00	-15.00 to +6.00	-15.00 to +6.00	
ADD POWER	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	
CYLINDER	to -7.00	to -7.00	to -7.00	to -7.00	
OTHER		Gray/Brown	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS - 11mm	16 250 67	T6 250 67	X6 250 67	PZ6 250 67	
INVISIBLE MARKINGS - 13mm	16 250 67	T6 250 67	X6 250 67	PZ6 250 67	
INVISIBLE MARKINGS - 15mm	16 250 67	T6 250 67	X6 250 67	PZ6 250 67	
INVISIBLE MARKINGS - 18mm	16 250 67	T6 250 67	X6 250 67	PZ6 250 67	
INVISIBLE MARKINGS - Variable	16 250 67	T6 250 67	X6 250 67	PZ6 250 67	

SuperLite™ 1.74	TREATMENT		COATING
	Clear	Transitions Signature7™	
POWER RANGE	-18.75 to +15.50	-18.75 to +15.50	
ADD POWER	+0.75 to +4.00	+0.75 to +4.00	
CYLINDER	to -9.00	to -9.00	
OTHER		Gray/Brown	
INVISIBLE MARKINGS - 11mm	74 250 67	T74 250 67	
INVISIBLE MARKINGS - 13mm	74 250 67	T74 250 67	
INVISIBLE MARKINGS - 15mm	74 250 67	T74 250 67	
INVISIBLE MARKINGS - 18mm	74 250 67	T74 250 67	
INVISIBLE MARKINGS - Variable	74 250 67	T74 250 67	

EVERYDAY> PROGRESSIVE

SHAMIR
Autograph II®

PRESCRIPTOR ORDER CODE:

AUTOGRAPH II 11
AUTOGRAPH II 13
AUTOGRAPH II 15
AUTOGRAPH II 18

PRODUCT FEATURES:

- ☒ FREEFORM PAL
☐ FREEFORM SV
☐ CONVENTIONAL PAL
☐ CONVENTIONAL SV
☒ EYEPOINT TECHNOLOGY®
☒ AS-WORN TECHNOLOGY™

NOTE:

- Recommended min fitting heights are 11mm, 13mm, 15mm and 18mm
- Prism to 10.00 Diopters

DLC™ TRIVEX™	TREATMENT					COATING
	Clear	Transitions Signature7™	XTRActive	Vantage	Polarized	
POWER RANGE	-13.25 to +6.00	-13.25 to +6.00	-13.25 to +6.00	-13.25 to +6.00	-13.25 to +6.00	
ADD POWER	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	
CYLINDER	to -6.50	to -6.50	to -6.50	to -6.50	to -6.50	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS - 11mm	TX 250 611	TTX 250 611	XX 250 611	VX 250 611	PTX 250 611	
INVISIBLE MARKINGS - 13mm	TX 250 613	TTX 250 613	XX 250 613	VX 250 613	PTX 250 613	
INVISIBLE MARKINGS - 15mm	TX 250 615	TTX 250 615	XX 250 615	VX 250 615	PTX 250 615	
INVISIBLE MARKINGS - 18mm	TX 250 618	TTX 250 618	XX 250 618	VX 250 618	PTX 250 618	
INVISIBLE MARKINGS - Variable	TX 250 61	TTX 250 61	XX 250 61	VX 250 61	PTX 250 61	

POLYPLUS™ / POLYCARBONATE	TREATMENT						COATING
	Clear	Transitions Signature7™	XTRActive	Vantage	Polarized	Drivewear	
POWER RANGE	-14.75 to +6.50	-14.75 to +6.50	-14.75 to +6.50	-14.75 to +6.50	-14.75 to +6.50	-14.75 to +6.50	
ADD POWER	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	
CYLINDER	to -7.00	to -7.00	to -7.00	to -7.00	to -7.00	to -7.00	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15		
INVISIBLE MARKINGS - 11mm	PC 250 611	PT 250 611	PX 250 611	PV 250 611	PPZ 250 611	PD 250 611	
INVISIBLE MARKINGS - 13mm	PC 250 613	PT 250 613	PX 250 613	PV 250 613	PPZ 250 613	PD 250 613	
INVISIBLE MARKINGS - 15mm	PC 250 615	PT 250 615	PX 250 615	PV 250 615	PPZ 250 615	PD 250 615	
INVISIBLE MARKINGS - 18mm	PC 250 618	PT 250 618	PX 250 618	PV 250 618	PPZ 250 618	PD 250 618	
INVISIBLE MARKINGS - Variable	PC 250 61	PT 250 61	PX 250 61	PV 250 61	PPZ 250 61	PD 250 61	

NXT	TREATMENT		COATING
	Photochromic	Polarized	
POWER RANGE	-1.50 to +7.00	-1.50 to +7.00	
ADD POWER	+0.75 to +4.00	+0.75 to +4.00	
CYLINDER	to -6.50	to -6.50	
OTHER	Green/Copper	Gray/Brown	
INVISIBLE MARKINGS - 11mm	NXT 250 611	NXTZ 250 611	
INVISIBLE MARKINGS - 13mm	NXT 250 613	NXTZ 250 613	
INVISIBLE MARKINGS - 15mm	NXT 250 615	NXTZ 250 615	
INVISIBLE MARKINGS - 18mm	NXT 250 618	NXTZ 250 618	
INVISIBLE MARKINGS - Variable	NXT 250 61	NXTZ 250 61	

Crizal AR Coating Hoya AR Coating Zeiss AR Coating iCoat AR Coating

- Denotes that material is available in a Shamir Freeform Blank™

All Shamir products are available with your lab's in-house AR Coatings. Please check with your lab to determine which specific manufacturer's AR coating is available with the Shamir lens and material that you are ordering, as a manufacturer may offer more than one AR coating.

EFFECTIVE DATE: SEP. 15, 2014



EVERYDAY> PROGRESSIVE

Shamir InTouch™

PRESCRIPTOR ORDER CODE:

INTOUCH 15

INTOUCH 18














PRODUCT FEATURES:










- ☒ FREEFORM PAL
- ☐ FREEFORM SV
- ☐ CONVENTIONAL PAL
- ☐ CONVENTIONAL SV










- ☒ EYEPOINT TECHNOLOGY™
- ☐ AS-WORN TECHNOLOGY™






NOTE:










- Recommended min fitting heights are 15mm and 18mm
- Prism to 10.00 Diopters






CR-39	TREATMENT						COATING
	Clear	Transitions Signature7™	XTRActive	Vantage	Polarized	Drivewear	
POWER RANGE	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	
ADD POWER	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	
CYLINDER	to -6.00	to -6.00	to -6.00	to -6.00	to -6.00	to -6.00	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15		
INVISIBLE MARKINGS - 15mm	 15 250 I15	 TR 250 I15	 X5 250 I15	 V5 250 I15	 PZ 250 I15	 DW 250 I15	
INVISIBLE MARKINGS - 18mm	 15 250 I18	 TR 250 I18	 X5 250 I18	 V5 250 I18	 PZ 250 I18	 DW 250 I18	














SuperLite™ 1.60	TREATMENT				COATING
	Clear	Transitions Signature7™	XTRActive	Polarized	
POWER RANGE	-15.00 to +6.00	-15.00 to +6.00	-15.00 to +6.00	-15.00 to +6.00	
ADD POWER	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	
CYLINDER	to -7.00	to -7.00	to -7.00	to -7.00	
OTHER		Gray/Brown	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS - 15mm	 16 250 I15	 T6 250 I15	 X6 250 I15	 PZ6 250 I15	
INVISIBLE MARKINGS - 18mm	 16 250 I18	 T6 250 I18	 X6 250 I18	 PZ6 250 I18	












SuperLite™ 1.67	TREATMENT				COATING
	Clear	Transitions Signature7™	XTRActive	Polarized	
POWER RANGE	-16.75 to +8.00	-16.75 to +8.00	-16.75 to +8.00	-16.75 to +8.00	
ADD POWER	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	
CYLINDER	to -8.00	to -8.00	to -8.00	to -8.00	
OTHER		Gray/Brown	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS - 15mm	 67 250 I15	 T67 250 I15	 X67 250 I15	 PZ67 250 I15	
INVISIBLE MARKINGS - 18mm	 67 250 I18	 T67 250 I18	 X67 250 I18	 PZ67 250 I18	





SuperLite™ 1.74	TREATMENT		COATING
	Clear	Transitions Signature7™	
POWER RANGE	-18.75 to +15.50	-18.75 to +15.50	
ADD POWER	+0.75 to +3.50	+0.75 to +3.50	
CYLINDER	to -9.00	to -9.00	
OTHER		Gray/Brown	
INVISIBLE MARKINGS - 15mm	 74 250 I15	 T74 250 I15	
INVISIBLE MARKINGS - 18mm	 74 250 I18	 T74 250 I18	

DLC™ TRIVEX™	TREATMENT				COATING
	Clear	Transitions Signature7™	XTRActive	Polarized	
POWER RANGE	-13.25 to +6.00	-13.25 to +6.00	-13.25 to +6.00	-13.25 to +6.00	
ADD POWER	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	
CYLINDER	to -6.50	to -6.50	to -6.50	to -6.50	
OTHER		Gray/Brown	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS - 15mm	 TX 250 I15	 TTX 250 I15	 XX 250 I15	 PTX 250 I15	
INVISIBLE MARKINGS - 18mm	 TX 250 I18	 TTX 250 I18	 XX 250 I18	 PTX 250 I18	

1.56	TREATMENT		COATING
	Clear	Polarized	
POWER RANGE	-14.00 to +6.00	-14.00 to +6.00	
ADD POWER	+0.75 to +3.50	+0.75 to +3.50	
CYLINDER	to -6.75	to -6.75	
OTHER			
INVISIBLE MARKINGS - 15mm	 56 250 I15	 56P 250 I15	
INVISIBLE MARKINGS - 18mm	 56 250 I18	 56P 250 I18	


POLYPLUS™/ POLYCARBONATE	TREATMENT						COATING
	Clear	Transitions Signature7™	XTRActive	Vantage	Polarized	Drivewear	
POWER RANGE	-14.75 to +6.00	-14.75 to +6.00	-14.75 to +6.00	-14.75 to +6.00	-14.75 to +6.00	-14.75 to +6.00	
ADD POWER	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	
CYLINDER	to -7.00	to -7.00	to -7.00	to -7.00	to -7.00	to -7.00	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15		
INVISIBLE MARKINGS - 15mm	 PC 250 I15	 PT 250 I15	 PX 250 I15	 PV 250 I15	 PPZ 250 I15	 PT 250 I15	
INVISIBLE MARKINGS - 18mm	 PC 250 I18	 PT 250 I18	 PX 250 I18	 PV 250 I18	 PPZ 250 I18	 PT 250 I18	

Mineral 1.52	TREATMENT					COATING
	White	Photochromic	Tints	Polarized	Polarized & Photochromic	
POWER RANGE	-13.25 to 6.00	-13.25 to 6.00	-13.25 to 6.00	-13.25 to 6.00	-13.25 to 6.00	
ADD POWER	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	
CYLINDER	to -6.50	to -6.50	to -6.50	to -6.50	to -6.50	
OTHER		Thin&Dark	*	**	Gray/Brown/Autumn Gold	
INVISIBLE MARKINGS - 15mm	 G15 250 I15	 G15 250 I15	 G15 250 I15	 G15 250 I15	 G15 250 I15	
INVISIBLE MARKINGS - 18mm	 G15 250 I18	 G15 250 I18	 G15 250 I18	 G15 250 I18	 G15 250 I18	

 Crizal AR Coating
  Hoya AR Coating
  Zeiss AR Coating
  iCoat AR Coating

All Shamir products are available with your lab's in-house AR Coatings. Please check with your lab to determine which specific manufacturer's AR coating is available with the Shamir lens and material that you are ordering, as a manufacturer may offer more than one AR coating.

EFFECTIVE DATE: SEP. 15, 2014

 - Denotes that material is available in a Shamir Freeform Blank™

* Tints- Autumn Gold, Canary Yellow, Rose 1/2, Gray 1/2/3, Tan 3, Green 1/2/3, G-15

** Polarized- Autumn Gold, Green, Gray 1/3, Brown 1/3, G15



SHAMIR

EVERYDAY> PROGRESSIVE

Shamir InTouch™





PREScriPTOR ORDER CODE:





INTOUCH 15 INTOUCH 18

- PRODUCT FEATURES:
- ☒ FREEFORM PAL
☐ FREEFORM SV
☐ CONVENTIONAL PAL
☐ CONVENTIONAL SV
- ☒ EYEPOINT TECHNOLOGY®
☐ AS-WORN TECHNOLOGY™

NOTE:


- Recommended min fitting heights are 15mm and 18mm
- Prism to 10.00 Diopters

Mineral 1.6	TREATMENT		COATING
	White		
	POWER RANGE	-15.00 to +6.00	
	ADD POWER	+0.75 to +3.50	
	CYLINDER	to -7.00	
	OTHER		
	INVISIBLE MARKINGS - 15mm	 250  G16 I15	
INVISIBLE MARKINGS - 18mm	 250  G16 I18		

 Crizal AR Coating  Hoya AR Coating  Zeiss AR Coating  iCoat AR Coating

All Shamir products are available with your lab's in-house AR Coatings. Please check with your lab to determine which specific manufacturer's AR coating is available with the Shamir lens and material that you are ordering, as a manufacturer may offer more than one AR coating.

EFFECTIVE DATE: SEP.15, 2014

 - Denotes that material is available in a Shamir Freeform Blank™

* Tints- Autumn Gold, Canary Yellow, Rose 1/2, Gray 1/2/3, Tan 3, Green 1/2/3, G-15
** Polarized- Autumn Gold, Green, Gray 1/3, Brown 1/3, G15





PREScriPTOR ORDER CODE:
SPECTRUM 14 SPECTRUM 18
SPECTRUM 16

- PRODUCT FEATURES:
- ☒ FREEFORM PAL
 - ☐ FREEFORM SV
 - ☐ CONVENTIONAL PAL
 - ☐ CONVENTIONAL SV
 - ☒ EYEPOINT TECHNOLOGY®
 - ☐ AS-WORN TECHNOLOGY™

NOTE:

- Recommended min fitting heights are 14mm, 16mm and 18mm
- Prism to 10.00 Diopters

CR-39	TREATMENT						COATING
	Clear	Transitions Signature7™	XTRActive	Vantage	Polarized	Drivewear	
POWER RANGE	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	
ADD POWER	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	
CYLINDER	to -6.00	to -6.00	to -6.00	to -6.00	to -6.00	to -6.00	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15		
INVISIBLE MARKINGS - 14mm	15 250 SP14	TR 250 SP14	X5 250 SP14	V5 250 SP14	PZ 250 SP14	DW 250 SP14	
INVISIBLE MARKINGS - 16mm	15 250 SP16	TR 250 SP16	X5 250 SP16	V5 250 SP16	PZ 250 SP16	DW 250 SP16	
INVISIBLE MARKINGS - 18mm	15 250 SP18	TR 250 SP18	X5 250 SP18	V5 250 SP18	PZ 250 SP18	DW 250 SP18	

SuperLite™ 1.60	TREATMENT				COATING
	Clear	Transitions Signature7™	XTRActive	Polarized	
POWER RANGE	-15.00 to +6.00	-15.00 to +6.00	-15.00 to +6.00	-15.00 to +6.00	
ADD POWER	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	
CYLINDER	to -7.00	to -7.00	to -7.00	to -7.00	
OTHER		Gray/Brown	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS - 14mm	16 250 SP14	T6 250 SP14	X6 250 SP14	PZ6 250 SP14	
INVISIBLE MARKINGS - 16mm	16 250 SP16	T6 250 SP16	X6 250 SP16	PZ6 250 SP16	
INVISIBLE MARKINGS - 18mm	16 250 SP18	T6 250 SP18	X6 250 SP18	PZ6 250 SP18	

SuperLite™ 1.67	TREATMENT				COATING
	Clear	Transitions Signature7™	XTRActive	Polarized	
POWER RANGE	-16.75 to +8.00	-16.75 to +8.00	-16.75 to +8.00	-16.75 to +8.00	
ADD POWER	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	
CYLINDER	to -8.00	to -8.00	to -8.00	to -8.00	
OTHER		Gray/Brown	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS - 14mm	67 250 SP14	T67 250 SP14	X67 250 SP14	PZ67 250 SP14	
INVISIBLE MARKINGS - 16mm	67 250 SP16	T67 250 SP16	X67 250 SP16	PZ67 250 SP16	
INVISIBLE MARKINGS - 18mm	67 250 SP18	T67 250 SP18	X67 250 SP18	PZ67 250 SP18	

SuperLite™ 1.74	TREATMENT		COATING
	Clear	Transitions Signature7™	
POWER RANGE	-18.75 to +14.75	-18.75 to +14.75	
ADD POWER	+0.75 to +3.50	+0.75 to +3.50	
CYLINDER	to -9.00	to -9.00	
OTHER			
INVISIBLE MARKINGS - 14mm	74 250 SP14	T74 250 SP14	
INVISIBLE MARKINGS - 16mm	74 250 SP16	T74 250 SP16	
INVISIBLE MARKINGS - 18mm	74 250 SP18	T74 250 SP18	

DLC™ TRIVEX™	TREATMENT					COATING
	Clear	Transitions Signature7™	XTRActive	Vantage	Polarized	
POWER RANGE	-13.25 to +6.00	-13.25 to +6.00	-13.25 to +6.00	-13.25 to +6.00	-13.25 to +6.00	
ADD POWER	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	
CYLINDER	to -6.50	to -6.50	to -6.50	to -6.50	to -6.50	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS - 14mm	TX 250 SP14	TTX 250 SP14	XX 250 SP14	VX 250 SP14	PTX 250 SP14	
INVISIBLE MARKINGS - 16mm	TX 250 SP16	TTX 250 SP16	XX 250 SP16	VX 250 SP16	PTX 250 SP16	
INVISIBLE MARKINGS - 18mm	TX 250 SP18	TTX 250 SP18	XX 250 SP18	VX 250 SP18	PTX 250 SP18	

POLYPLUS™/ POLYCARBONATE	TREATMENT						COATING
	Clear	Transitions Signature7™	XTRActive	Vantage	Polarized	Drivewear	
POWER RANGE	-14.75 to +6.00	-14.75 to +6.00	-14.75 to +6.00	-14.75 to +6.00	-14.75 to +6.00	-14.75 to +6.00	
ADD POWER	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	
CYLINDER	to -7.00	to -7.00	to -7.00	to -7.00	to -7.00	to 7.00	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15		
INVISIBLE MARKINGS - 14mm	PC 250 SP14	PT 250 SP14	PX 250 SP14	PV 250 SP14	PPZ 250 SP14	PD 250 SP14	
INVISIBLE MARKINGS - 16mm	PC 250 SP16	PT 250 SP16	PX 250 SP16	PV 250 SP16	PPZ 250 SP16	PD 250 SP16	
INVISIBLE MARKINGS - 18mm	PC 250 SP18	PT 250 SP18	PX 250 SP18	PV 250 SP18	PPZ 250 SP18	PD 250 SP18	





PREScriPTOR ORDER CODE:











SPECTRUM 14 SPECTRUM 18
SPECTRUM 16







PRODUCT FEATURES:




















- ☒ FREEFORM PAL
☐ FREEFORM SV
☐ CONVENTIONAL PAL
☐ CONVENTIONAL SV
- ☒ EYEPOINT TECHNOLOGY®
☐ AS-WORN TECHNOLOGY™

NOTE:

- Recommended min fitting heights are 14mm, 16mm and 18mm
- Prism to 10.00 Diopters

1.56	TREATMENT		COATING
	Clear	Polarized	
POWER RANGE	-14.00 to +6.00	-14.00 to +6.00	   
ADD POWER	+0.75 to +3.50	+0.75 to +3.50	
CYLINDER	to -6.75	to -6.75	
OTHER			
INVISIBLE MARKINGS - 14mm	 56 250 SP14	 56P 250 SP14	
INVISIBLE MARKINGS - 16mm	 56 250 SP16	 56P 250 SP16	
INVISIBLE MARKINGS - 18mm	 56 250 SP18	 56P 250 SP18	

Mineral 1.60	TREATMENT		COATING
	White		
	POWER RANGE	-15.00 to +6.00	
	ADD POWER	+0.75 to +3.50	
	CYLINDER	to -7.00	
	OTHER		
	INVISIBLE MARKINGS - 14mm	 250  G16 SP14	
	INVISIBLE MARKINGS - 16mm	 250  G16 SP16	
INVISIBLE MARKINGS - 18mm	 250  G16 SP18		

Mineral 1.52	TREATMENT					COATING
	White	Photochromic	Tints	Polarized	Polarized & Photochromic	
POWER RANGE	-13.25 to +6.00	-13.25 to +6.00	-13.25 to +6.00	-13.25 to +6.00	-13.25 to +6.00	   
ADD POWER	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	
CYLINDER	to -6.50	to -6.50	to -6.50	to -6.50	to -6.50	
OTHER		Thin&Dark	*	**	Gray/Brown/ Autumn Gold	
INVISIBLE MARKINGS - 14mm	 250	 G15 SP14	 G15 SP14	 GP15 SP14	 G15 SP14	
INVISIBLE MARKINGS - 16mm	 250	 G15 SP16	 G15 SP16	 GP15 SP16	 G15 SP16	
INVISIBLE MARKINGS - 18mm	 250	 G15 SP18	 G15 SP18	 GP15 SP18	 G15 SP18	

PRESCRIPTOR ORDER CODE:

ELEMENT 16

ELEMENT 19











PRODUCT FEATURES:











- ☒ FREEFORM PAL
☐ FREEFORM SV
☐ CONVENTIONAL PAL
☐ CONVENTIONAL SV









- ☒ EYEPOINT TECHNOLOGY®
☐ AS-WORN TECHNOLOGY™













NOTE:











- Recommended min fitting heights are 16mm and 19mm
- Prism to 10.00 Diopters

CR-39	TREATMENT			COATING
	Clear	Transitions Signature7™	Polarized	
POWER RANGE	-7.00 to +4.00	-7.00 to +4.00	-7.00 to +4.00	   
ADD POWER	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	
CYLINDER	to -6.00	to -6.00	to -6.00	
OTHER		Gray/Brown	Gray/Brown/G-15	
INVISIBLE MARKINGS - 16mm	 15 ES 250	 TR ES 250	 PZ ES 250	
INVISIBLE MARKINGS - 19mm	 15 E 250	 TR E 250	 PZ E 250	

DLC™ TRIVEX™	TREATMENT			COATING
	Clear	Transitions Signature7™	XTRActive	
POWER RANGE	-13.25 to +6.00	-13.25 to +6.00	-13.25 to +6.00	   
ADD POWER	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	
CYLINDER	to -6.50	to -6.50	to -6.50	
OTHER		Gray/Brown	Gray	
INVISIBLE MARKINGS - 16mm	 TX ES 250	 TTX ES 250	 XX ES 250	
INVISIBLE MARKINGS - 19mm	 TX E 250	 TTX E 250	 XX E 250	

SuperLite™ 1.60	TREATMENT		COATING
	Clear	Transitions Signature7™	
POWER RANGE	-10.00 to +6.00	-10.00 to +6.00	   
ADD POWER	+0.75 to +3.50	+0.75 to +3.50	
CYLINDER	to -7.00	to -7.00	
OTHER		Gray/Brown	
INVISIBLE MARKINGS - 16mm	 16 ES 250	 T6 ES 250	
INVISIBLE MARKINGS - 19mm	 16 E 250	 T6 E 250	

POLYPLUS™/ POLYCARBONATE	TREATMENT				COATING
	Clear	Transitions Signature7™	XTRActive	Polarized	
POWER RANGE	-10.00 to +6.00	-10.00 to +6.00	-10.00 to +6.00	-10.00 to +6.00	   
ADD POWER	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	
CYLINDER	to -7.00	to -7.00	to -7.00	to -7.00	
OTHER		Gray/Brown	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS - 16mm	 PC ES 250	 PT ES 250	 PX ES 250	 PPZ ES 250	
INVISIBLE MARKINGS - 19mm	 PC E 250	 PT E 250	 PX E 250	 PPZ E 250	














SuperLite™ 1.67	TREATMENT			COATING
	Clear	Transitions Signature7™	XTRActive	
POWER RANGE	-12.50 to +8.00	-12.50 to +8.00	-12.50 to +8.00	   
ADD POWER	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	
CYLINDER	to -8.00	to -8.00	to -8.00	
OTHER		Gray/Brown	Gray	
INVISIBLE MARKINGS - 16mm	 67 ES 250	 T67 ES 250	 X67 ES 250	
INVISIBLE MARKINGS - 19mm	 67 E 250	 T67 E 250	 X67 E 250	













PREScriPTOR ORDER CODE:
FIRSTPAL














- PRODUCT FEATURES:
- ☒ FREEFORM PAL
 - ☐ FREEFORM SV
 - ☐ CONVENTIONAL PAL
 - ☐ CONVENTIONAL SV
 - ☒ EYEPOINT TECHNOLOGY®
 - ☐ AS-WORN TECHNOLOGY™





NOTE:





- Recommended min fitting height is 15mm
- Prism to 10.00 Diopters







CR-39	TREATMENT					COATING
	 Clear	 Transitions Signature7™	XTRActive	 Vantage	 Polarized	
POWER RANGE	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	   
ADD POWER	+0.75 to +1.50	+0.75 to +1.50	+0.75 to +1.50	+0.75 to +1.50	+0.75 to +1.50	
CYLINDER	to -6.00	to -6.00	to -6.00	to -6.00	to -6.00	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS	 150 FP	 150 TR FP	 150 X5 FP	 150 V5 FP	 150 PZ FP	

SuperLite™ 1.60	TREATMENT				COATING
	 Clear	 Transitions Signature7™	 XTRActive	 Polarized	
POWER RANGE	-15.00 to +6.00	-15.00 to +6.00	-15.00 to +6.00	-15.00 to +6.00	   
ADD POWER	+0.75 to +1.50	+0.75 to +1.50	+0.75 to +1.50	+0.75 to +1.50	
CYLINDER	to -7.00	to -7.00	to -7.00	to -7.00	
OTHER		Gray/Brown	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS	 150 16 FP	 150 T6 FP	 150 X6 FP	 150 PZ6 FP	

POLYPLUS™/ POLYCARBONATE	TREATMENT					COATING
	 Clear	 Transitions Signature7™	XTRActive	 Vantage	 Polarized	
POWER RANGE	-14.75 to +6.00	-14.75 to +6.00	-14.75 to +6.00	-14.75 to +6.00	-14.75 to +6.00	   
ADD POWER	+0.75 to +1.50	+0.75 to +1.50	+0.75 to +1.50	+0.75 to +1.50	+0.75 to +1.50	
CYLINDER	to -7.00	to -7.00	to -7.00	to -7.00	to -7.00	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS	 150 PC FP	 150 PT FP	 150 PX FP	 150 PV FP	 150 PPZ FP	

CR-39	TREATMENT				
	Clear	Transitions®	XTRActive	Vantage	Polarized
POWER RANGE	-8.50 to +6.75	-8.50 to +6.75	--	--	--
ADD POWER	+0.75 to +3.50	+0.75 to +3.50	--	--	--
NOMINAL BASE CURVE	2, 3, 4.5, 5.5, 6.5, 7.5	2, 3, 4.5, 5.5, 6.5, 7.5	--	--	--
OTHER/COATING	 Hard Coated/Uncoated	 Gray/Brown	--	--	--
INVISIBLE MARKINGS	 250	 250	--	--	--

SuperLite™ 1.60	TREATMENT				
	Clear	Transitions®	XTRActive	Vantage	Polarized
POWER RANGE	-12.00 to +8.00	-12.00 to +8.00	--	--	--
ADD POWER	+1.00 to +3.50	+1.00 to +3.50	--	--	--
NOMINAL BASE CURVE	1.75, 3, 4.5, 5.5, 6.5, 7.5	1.75, 3, 4.5, 5.5, 6.5, 7.5	--	--	--
OTHER/COATING	 Hard Coated/Uncoated	 Gray	--	--	--
INVISIBLE MARKINGS	 H 250	 H 250	--	--	--

SuperLite™ 1.67	TREATMENT				
	Clear	Transitions®	XTRActive	Vantage	Polarized
POWER RANGE	-12.50 to +8.00	-12.50 to +8.00	--	--	-12.50 to +8.00
ADD POWER	+1.00 to +3.50	+1.00 to +3.50	--	--	+1.00 to +3.50
NOMINAL BASE CURVE	1.5, 3, 4.5, 5.5, 6.5, 7.5	1.5, 3, 4.5, 5.5, 6.5, 7.5	--	--	1.5, 3, 4.5, 5.5, 6.5, 7.5
OTHER/COATING	 Hard Coated/Uncoated	 Gray	--	--	 Gray
INVISIBLE MARKINGS	 67 250	 67 250	--	--	 167 250

PRODUCT FEATURES:






- ☐ FREEFORM PAL
- ☐ FREEFORM SV
- ☒ CONVENTIONAL PAL
- ☐ CONVENTIONAL SV
- ☒ EYEPOINT TECHNOLOGY®
- ☐ AS-WORN TECHNOLOGY®

NOTE:

- Recommended min fitting height is 19mm
- Cylinder to -6.00
- Prism to 4.00 Diopters

SuperLite™ 1.74	TREATMENT				
	Clear	Transitions®	XTRActive	Vantage	Polarized
POWER RANGE	--	--	--	--	--
ADD POWER	--	--	--	--	--
NOMINAL BASE CURVE	--	--	--	--	--
OTHER/COATING	--	--	--	--	--
INVISIBLE MARKINGS	--	--	--	--	--

DLC™ TRIVEX™	TREATMENT			
	Clear	Transitions®	XTRActive	Polarized
POWER RANGE	--	--	--	--
ADD POWER	--	--	--	--
NOMINAL BASE CURVE	--	--	--	--
OTHER/COATING	--	--	--	--
INVISIBLE MARKINGS	--	--	--	--

POLYPLUS™/ POLYCARBONATE	TREATMENT				
	Clear	Transitions®	XTRActive	Vantage	Polarized
POWER RANGE	-12.00 to +7.50	-12.00 to +7.50	--	--	--
ADD POWER	+1.00 to +3.50	+1.00 to +3.50	--	--	--
NOMINAL BASE CURVE	1.75, 3, 4.5, 5.5, 6.5, 7.5	1.75, 3, 4.5, 5.5, 6.5, 7.5	--	--	--
OTHER/COATING	Hard Coated:  Non Primer: 	 Gray	--	--	--
INVISIBLE MARKINGS	 P 250	 P 250	--	--	--



CR-39	TREATMENT				
	Clear	Transitions	XTRActive	Vantage	Polarized
POWER RANGE	-8.50 to +6.00	-8.50 to +6.00	--	--	--
ADD POWER	+0.75 to +3.50	+0.75 to +3.50	--	--	--
BASE CURVE	1.5, 3, 5, 6.5, 7.5	1.5, 3, 5, 6.5, 7.5	--	--	--
OTHER/COATING	Hard Coated/Uncoated	Gray	--	--	--
INVISIBLE MARKINGS	250	250	--	--	--

SuperLite™ 1.60	TREATMENT				
	Clear	Transitions	XTRActive	Vantage	Polarized
POWER RANGE	-11.25 to +7.00	-11.25 to +7.00	--	--	--
ADD POWER	+0.75 to +3.50	+1.00 to +3.50	--	--	--
BASE CURVE	1.5, 3, 5, 6.5, 7.5	1.5, 3, 5, 6.5, 7.5	--	--	--
OTHER/COATING	Hard Coated/Uncoated	Gray	--	--	--
INVISIBLE MARKINGS	H 250	H 250	--	--	--

SuperLite™ 1.67	TREATMENT				
	Clear	Transitions	XTRActive	Vantage	Polarized
POWER RANGE	-12.50 to +6.00	-12.50 to +6.00	--	--	--
ADD POWER	+1.00 to +3.00	+1.00 to +3.00	--	--	--
BASE CURVE	1.5, 3, 4.5, 5.75, 7	1.5, 3, 4.5, 5.75, 7	--	--	--
OTHER/COATING	Hard Coated/Uncoated	Gray	--	--	--
INVISIBLE MARKINGS	67 250	67 250	--	--	--

PRODUCT FEATURES:

- ☐ FREEFORM PAL
- ☐ FREEFORM SV
- ☒ CONVENTIONAL PAL
- ☐ CONVENTIONAL SV
- ☒ EYEPOINT TECHNOLOGY®
- ☐ AS-WORN TECHNOLOGY®

NOTE:







- Recommended min fitting height is 16mm
- Cylinder to -6.00
- Prism to 4.00 Diopters





SuperLite™ 1.74	TREATMENT				
	Clear				
POWER RANGE	--				
ADD POWER	--				
NOMINAL BASE CURVE	--				
OTHER/COATING	--				
INVISIBLE MARKINGS	--				



DLC™ TRIVEX™	TREATMENT			
	Clear	Transitions	XTRActive	Polarized
POWER RANGE	--	--	--	--
ADD POWER	--	--	--	--
NOMINAL BASE CURVE	--	--	--	--
OTHER/COATING	--	--	--	--
INVISIBLE MARKINGS	--	--	--	--

POLYPLUS™/ POLYCARBONATE	TREATMENT				
	Clear	Transitions	XTRActive	Vantage	Polarized
POWER RANGE	-8.50 to +6.00	-8.50 to +6.00	--	--	--
ADD POWER	+1.00 to +3.00	+1.00 to +3.00	--	--	--
BASE CURVE	1.5, 3, 5, 7.5	1.5, 3, 5, 7.5	--	--	--
OTHER/COATING	Hard Coated: Non Primer:	Gray	--	--	--
INVISIBLE MARKINGS	P 250	P 250	--	--	--

Glass	TREATMENT	
	Clear	Photo
POWER RANGE	-5.25 to +6.00	-5.25 to +6.00
ADD POWER	+0.75 to +3.50	+0.75 to +3.50
BASE CURVE	3, 5, 7.5	3, 5, 7.5
OTHER/COATING		
INVISIBLE MARKINGS	250	TD 250

CR-39	TREATMENT				
	Clear	Transitions®	XTRActive	Vantage	Polarized
POWER RANGE	-8.50 to +5.75	-8.50 to +5.75	--	--	-5.75 to +6.00
ADD POWER	+0.75 to +3.50	+0.75 to +3.50	--	--	+1.00 to +3.00
BASE CURVE	1.5, 3, 5, 6.5, 7.5	1.5, 3, 5, 6.5, 7.5	--	--	3, 5, 7.5
OTHER/COATING	 Hard Coated/Uncoated	 Gray	--	--	 Gray/Brown
INVISIBLE MARKINGS	 250	 250	--	--	 PZ 250

SuperLite™ 1.60	TREATMENT				
	Clear	Transitions®	XTRActive	Vantage	Polarized
POWER RANGE	-11.25 to +7.00	-11.25 to +7.00	--	--	--
ADD POWER	+0.75 to +3.50	+1.00 to +3.50	--	--	--
BASE CURVE	1.5, 3, 5, 6.5, 7.5	1.5, 3, 5, 7.5	--	--	--
OTHER/COATING	 Hard Coated/Uncoated	 Gray	--	--	--
INVISIBLE MARKINGS	 H 250	 H 250	--	--	--

SuperLite™ 1.67	TREATMENT				
	Clear	Transitions®	XTRActive	Vantage	Polarized
POWER RANGE	-12.50 to +7.50	--	--	--	--
ADD POWER	+1.00 to +3.50	--	--	--	--
BASE CURVE	1.5, 3, 5, 7.5	--	--	--	--
OTHER/COATING	 Hard Coated/Uncoated	--	--	--	--
INVISIBLE MARKINGS	 H67 250	--	--	--	--



PRODUCT FEATURES:



- ☐ FREEFORM PAL
- ☐ FREEFORM SV
- ☒ CONVENTIONAL PAL
- ☐ CONVENTIONAL SV
- ☒ EYEPOINT TECHNOLOGY®
- ☐ AS-WORN TECHNOLOGY®







NOTE:

- Recommended min fitting height is 16mm
- Cylinder to -6.00
- Prism to 4.00 Diopters

SuperLite™ 1.74	TREATMENT				
	Clear	Transitions®	XTRActive	Vantage	Polarized
POWER RANGE	--	--	--	--	--
ADD POWER	--	--	--	--	--
BASE CURVE	--	--	--	--	--
OTHER/COATING	--	--	--	--	--
INVISIBLE MARKINGS	--	--	--	--	--

DLC™ TRIVEX™	TREATMENT			
	Clear	Transitions®	XTRActive	Polarized
POWER RANGE	-5.75 to +6.00	--	--	--
ADD POWER	+1.00 to +3.00	--	--	--
BASE CURVE	3, 5, 7.5	--	--	--
OTHER/COATING	 Hard Coated/Uncoated	--	--	--
INVISIBLE MARKINGS	 TX 250	--	--	--

POLYPLUS™/ POLYCARBONATE	TREATMENT				
	Clear	Transitions®	XTRActive	Vantage	Polarized
POWER RANGE	-12.00 to +7.50	--	--	--	--
ADD POWER	+1.00 to +3.00	--	--	--	--
BASE CURVE	1.5, 3, 5, 7	--	--	--	--
OTHER/COATING	 PolyPlus™** - HC	--	--	--	--
INVISIBLE MARKINGS	 P 250	--	--	--	--

Glass	TREATMENT	
	Clear	Photo
POWER RANGE	-9.00 to +6.00	-9.00 to +6.00
ADD POWER	+0.75 to +3.50	+0.75 to +3.50
BASE CURVE	1.5, 3, 5, 7	1.5, 3, 5, 7
OTHER/COATING	 Clear	 Gray Extra
INVISIBLE MARKINGS	 16  250	 16  250

EVERYDAY> BIFOCAL

Shamir Duo™

PRESCRIBTOR ORDER CODE:














DUO



















PRODUCT FEATURES:
















- | | | |
|--|---|---|
| <input type="checkbox"/> FREEFORM PAL | <input type="checkbox"/> EYEPOINT TECHNOLOGY III® | <input type="checkbox"/> EYEPOINT TECHNOLOGY® |
| <input type="checkbox"/> FREEFORM OCCUPATION | <input type="checkbox"/> AS-WORN QUADRO™ | <input type="checkbox"/> AS-WORN TECHNOLOGY™ |
| <input type="checkbox"/> FREEFORM SV | <input type="checkbox"/> INTELLICORRIDOR™ | |
| <input checked="" type="checkbox"/> FREEFORM BIFOCAL | <input type="checkbox"/> NATURAL POSTURE™ | |
| <input type="checkbox"/> CONVENTIONAL PAL | | |
| <input type="checkbox"/> CONVENTIONAL OCCUPATION | | |
| <input type="checkbox"/> CONVENTIONAL SV | | |

NOTE:
















- Recommended min fitting height is 15mm
- Prism to 10.00 Diopters
















CR-39	TREATMENT						COATING
	 Clear	 Transitions Signature7™	XTRActive	 Vantage	 Polarized	Drivewear	
POWER RANGE	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	  
ADD POWER	0.75 to +3.50	0.75 to +3.50	0.75 to +3.50	0.75 to +3.50	0.75 to +3.50	0.75 to +3.50	
CYLINDER	to -6.00	to -6.00	to -6.00	to -6.00	to -6.00	to -6.00	
OTHER	Gray/Brown		Gray	Gray	Gray/Brown/G-15		
INVISIBLE MARKING	 15 250	 250	 250	 250	 250	 250	










DLC™/TRIVEX™	TREATMENT					COATING
	 Clear	 Transitions Signature7™	 XTRActive	 Vantage	 Polarized	
POWER RANGE	-13.25 to +6.00	-13.25 to +6.00	-13.25 to +6.00	-13.25 to +6.00	-13.25 to +6.00	  
ADD POWER	0.75 to +4.00	0.75 to +4.00	0.75 to +4.00	0.75 to +4.00	0.75 to +4.00	
CYLINDER	to -6.50	to -6.50	to -6.50	to -6.50	to -6.50	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15	
INVISIBLE MARKING	 TX 250 	 TTX 250 	 XX 250 	 VX 250 	 PTX 250 	

SuperLite™ 1.60	TREATMENT				COATING
	 Clear	 Transitions Signature7™	 XTRActive	 Polarized	
POWER RANGE	-15.00 to +6.00	-15.00 to +6.00	-15.00 to +6.00	-15.00 to +6.00	  
ADD POWER	0.75 to +3.50	0.75 to +3.50	0.75 to +3.50	0.75 to +3.50	
CYLINDER	to -7.00	to -7.00	to -7.00	to -7.00	
OTHER	Gray/Brown		Gray	Gray/Brown/G-15	
INVISIBLE MARKING	 16 350 	 T6 250 	 X6 250 	 PZ6 250 	

1.56	TREATMENT		COATING
	Clear	Polarized	
POWER RANGE	-14.00 to +6.00	-14.00 to +6.00	
ADD POWER	0.75 to +4.00	0.75 to +4.00	
CYLINDER	to -6.75	to -6.75	
OTHER	BluTech Indoor	BluTech Outdoor	
INVISIBLE MARKING	250	250	

SuperLite™ 1.67	TREATMENT				COATING
	 Clear	 Transitions Signature7™	 XTRActive	 Polarized	
POWER RANGE	-16.75 to +8.00	-16.75 to +8.00	-16.75 to +8.00	-16.75 to +8.00	
ADD RANGE	0.75 to +4.00	0.75 to +4.00	0.75 to +4.00	0.75 to +4.00	
CYLINDER	to -8.00	to -8.00	to -8.00	to -8.00	
OTHER	Gray/Brown		Gray	Gray/Brown/G-15	
INVISIBLE MARKING	 67 250 	 T67 250 	 X67 250 	 PZ67 250 	

POLYPLUS™/ POLYCARBONATE	TREATMENT						COATING
	 Clear	 Transitions Signature7™	 XTRActive	 Vantage	 Polarized	 Drivewear	
POWER RANGE	-14.75 to +6.00	-14.75 to +6.00	-14.75 to +6.00	-14.75 to +6.00	-14.75 to +6.00	-14.75 to +6.00	  
ADD POWER	0.75 to +4.00	0.75 to +4.00	0.75 to +4.00	0.75 to +4.00	0.75 to +4.00	0.75 to +4.00	
CYLINDER	to -7.00	to -7.00	to -7.00	to -7.00	to -7.00	to -7.00	
OTHER	Gray/Brown		Gray	Gray	Gray/Brown/G-15		
INVISIBLE MARKING	 PC	 PT	 PX	 PVT	 PPZ	 PD	

SuperLite™ 1.74	TREATMENT		COATING
	 Clear	 Transitions Signature7™	
POWER RANGE	-18.75 to +12.00	-18.75 to +12.00	
ADD RANGE	0.75 to +3.50	0.75 to +3.50	
CYLINDER	to -9.00	to -9.00	
OTHER	Gray/Brown		
INVISIBLE MARKING	 74  250  JA	 74  250  JA	

NXT	TREATMENT		COATING
	Photochromic	Polarized	
POWER RANGE	-1.75 to +6.00	-1.75 to +6.00	
ADD POWER	0.75 to +4.00	0.75 to +4.00	
CYLINDER	to -6.50	to -6.50	
OTHER	Green/Copper	Gray/Brown	
INVISIBLE MARKINGS	250	250	

Crizal AR Coating
 Hoya AR Coating
 Zeiss AR Coating
 iCoat AR Coating

All Shamir products are available with your lab's in-house AR Coatings. Please check with your lab to determine which specific manufacturer's AR coating is available with the Shamir lens and material that you are ordering, as a manufacturer may offer more than one AR coating.

EFFECTIVE DATE: SEP. 15, 2014

- Denotes that material is available in a Shamir Freeform Blank™

Shamir Autograph III® - SV

PRESCRIBTOR ORDER CODE:

AUTOGRAPH III SV

PRODUCT FEATURES:

- ☐ FREEFORM PAL
☐ FREEFORM OCCUPATION
☒ FREEFORM SV
☐ FREEFORM BIFOCAL
☐ CONVENTIONAL PAL
☐ CONVENTIONAL OCCUPATION
☐ CONVENTIONAL SV
- ☒ EYEPOINT TECHNOLOGY III®
☒ AS-WORN QUADRO™
☐ INTELLICORRIDOR™
☐ NATURAL POSTURE™
- ☐ EYEPOINT TECHNOLOGY®
☐ AS-WORN TECHNOLOGY™

NOTE:

• Prism to 10.00 Diopters

CR-39	TREATMENT						COATING
	Clear	Transitions Signature7™	XTRActive	Vantage	Polarized	Drivewear	
POWER RANGE	-12.00 to +7.25	-12.00 to +7.25	-12.00 to +7.25	-12.00 to +7.25	-12.00 to +7.25	-12.00 to +7.25	
ADD POWER							
CYLINDER	to -6.00	to -6.00	to -6.00	to -6.00	to -6.00	to -6.00	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15		
INVISIBLE MARKINGS	SV 15 III	SV TR III	SV X5 III	SV V5 III	SV PZ III	SV DW III	

SuperLite™ 1.74	TREATMENT		COATING
	Clear	Transitions Signature7™	
POWER RANGE	-20.00 to +13.50	-20.00 to +13.50	
ADD POWER			
CYLINDER	to -9.00	to -9.00	
OTHER		Gray/Brown	
INVISIBLE MARKINGS	SV 74 III	SV T74 III	

SuperLite™ 1.60	TREATMENT				COATING
	Clear	Transitions Signature7™	XTRActive	Polarized	
POWER RANGE	-14.00 to +8.75	-14.00 to +8.75	-14.00 to +8.75	-14.00 to +8.75	
ADD POWER					
CYLINDER	to -7.00	to -7.00	to -7.00	to -7.00	
OTHER		Gray/Brown	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS	SV 16 III	SV T6 III	SV X6 III	SV PZ6 III	

1.56	TREATMENT		COATING
	Clear	Polarized	
POWER RANGE	-14.00 to +8.50	-14.00 to +8.50	
ADD POWER			
CYLINDER	to -6.75	to -6.75	
OTHER	BluTech Indoor	BluTech Outdoor	
INVISIBLE MARKINGS	SV 56 III	SV 56P III	

SuperLite™ 1.67	TREATMENT				COATING
	Clear	Transitions Signature7™	XTRActive	Polarized	
POWER RANGE	-16.00 to +16.00	-16.00 to +16.00	-16.00 to +16.00	-16.00 to +16.00	
ADD POWER					
CYLINDER	to -8.00	to -8.00	to -8.00	to -8.00	
OTHER		Gray/Brown	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS	SV 67 III	SV T67 III	SV X67 III	SV PZ67 III	

POLYPLUS™/ POLYCARBONATE	TREATMENT						COATING
	Clear	Transitions Signature7™	XTRActive	Vantage	Polarized	Drivewear	
POWER RANGE	-14.00 to +8.75	-14.00 to +8.75	-14.00 to +8.75	-14.00 to +8.75	-14.00 to +8.75	-14.00 to +8.75	
ADD POWER							
CYLINDER	to -7.00	to -7.00	to -7.00	to -7.00	to -7.00	to -7.00	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15		
INVISIBLE MARKINGS	SV PC III	SV PT III	SV PX III	SV PV III	SV PPZ III	SV PD III	

DLC™ TRIVEX™	TREATMENT					COATING
	Clear	Transitions Signature7™	XTRActive	Vantage	Polarized	
POWER RANGE	-14.00 to +8.00	-14.00 to +8.00	-14.00 to +8.00	-14.00 to +8.00	-14.00 to +8.00	
ADD POWER						
CYLINDER	to -6.50	to -6.50	to -6.50	to -6.50	to -6.50	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G15	
INVISIBLE MARKINGS	SV TX III	SV TTX III	SV XX III	SV VX III	SV PTX III	

NXT	TREATMENT		COATING
	Photochromic	Polarized	
POWER RANGE	+1.00 to +8.00	+1.00 to +8.00	
ADD POWER			
CYLINDER	to -6.50	to -6.50	
OTHER	Green/Copper	Gray/Brown	
INVISIBLE MARKINGS	SV TX III	SV TX III	

Crizal AR Coating
 Hoya AR Coating
 Zeiss AR Coating
 iCoat AR Coating

All Shamir products are available with your lab's in-house AR Coatings. Please check with your lab to determine which specific manufacturer's AR coating is available with the Shamir lens and material that you are ordering, as a manufacturer may offer more than one AR coating.

EFFECTIVE DATE: FEB. 2, 2015

- Denotes that material is available in a Shamir Freeform Blank™

EVERYDAY> SINGLE VISION

SHAMIR

Autograph II SINGLE VISION™

PRESCRIPTOR ORDER CODE:

AUTOGRAPH II SV AW

PRODUCT FEATURES:

- ☐ FREEFORM PAL
- ☒ FREEFORM SV
- ☐ CONVENTIONAL PAL
- ☐ CONVENTIONAL SV
- ☒ EYEPOINT TECHNOLOGY®
- ☒ AS-WORN TECHNOLOGY™

NOTE:

• Prism to 10.00 Diopters

CR-39	TREATMENT						COATING
	Clear	Transitions Signature7™	XTRActive	Vantage	Polarized	Drivewear	
POWER RANGE	-14.00 to +8.00	-14.00 to +8.00	-14.00 to +8.00	-14.00 to +8.00	-14.00 to +8.00	-14.00 to +8.00	
ADD POWER							
CYLINDER	to -6.00	to -6.00	to -6.00	to -6.00	to -6.00	to -6.00	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15		
INVISIBLE MARKINGS	15 FF SV	TR FF SV	X5 FF SV	V5 FF SV	PZ FF SV	DW FF SV	

SuperLite™ 1.60	TREATMENT				COATING
	Clear	Transitions Signature7™	XTRActive	Polarized	
POWER RANGE	-14.00 to +8.00	-14.00 to +8.00	-14.00 to +8.00	-14.00 to +8.00	
ADD POWER					
CYLINDER	to -7.00	to -7.00	to -7.00	to -7.00	
OTHER		Gray/Brown	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS	16 FF SV	T6 FF SV	X6 FF SV	PZ6 FF SV	

SuperLite™ 1.67	TREATMENT				COATING
	Clear	Transitions Signature7™	XTRActive	Polarized	
POWER RANGE	-16.00 to +16.00	-16.00 to +16.00	-16.00 to +16.00	-16.00 to +16.00	
ADD POWER					
CYLINDER	to -8.00	to -8.00	to -8.00	to -8.00	
OTHER		Gray/Brown	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS	67 FF SV	T67 FF SV	X67 FF SV	PZ67 FF SV	

SuperLite™ 1.74	TREATMENT		COATING
	Clear	Transitions Signature7™	
POWER RANGE	-20.00 to +15.50	-20.00 to +15.50	
ADD POWER			
CYLINDER	to -9.00	to -9.00	
OTHER		Gray/Brown	
INVISIBLE MARKINGS	74 FF SV	T74 FF SV	

DLC™ TRIVEX™	TREATMENT					COATING
	Clear	Transitions Signature7™	XTRActive	Vantage	Polarized	
POWER RANGE	-14.00 to +7.50	-14.00 to +7.50	-14.00 to +7.50	-14.00 to +7.50	-14.00 to +7.50	
ADD POWER						
CYLINDER	to -6.50	to -6.50	to -6.50	to -6.50	to -6.50	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS	TX FF SV	TTX FF SV	TTX FF SV	VX FF SV	PTX FF SV	

1.56	TREATMENT		COATING
	Clear	Polarized	
POWER RANGE	-14.00 to +8.50	-14.00 to +8.50	
ADD POWER			
CYLINDER	to -6.50	to -6.50	
OTHER			
INVISIBLE MARKINGS	56 FF SV	56P FF SV	

POLYPLUS™/ POLYCARBONATE	TREATMENT						COATING
	Clear	Transitions Signature7™	XTRActive	Vantage	Polarized	Drivewear	
POWER RANGE	-14.00 to +8.75	-14.00 to +8.75	-14.00 to +8.75	-14.00 to +8.75	-14.00 to +8.75	-14.00 to +8.75	
ADD POWER							
CYLINDER	to -7.00	to -7.00	to -7.00	to -7.00	to -7.00	to -7.00	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15		
INVISIBLE MARKINGS	PC FF SV	PT FF SV	PX FF SV	PV FF SV	PPZ FF SV	PD FF SV	

Crizal AR Coating
 Hoya AR Coating
 Zeiss AR Coating
 iCoat AR Coating

All Shamir products are available with your lab's in-house AR Coatings. Please check with your lab to determine which specific manufacturer's AR coating is available with the Shamir lens and material that you are ordering, as a manufacturer may offer more than one AR coating.

EFFECTIVE DATE: SEP. 15, 2014

- Denotes that material is available in a Shamir Freeform Blank™

SHAMIR

EVERYDAY> SINGLE VISION



PRESCRIPTOR ORDER CODE:

SPECTRUM SV

PRODUCT FEATURES:

- ☐ FREEFORM PAL
- ☒ FREEFORM SV
- ☐ CONVENTIONAL PAL
- ☐ CONVENTIONAL SV
- ☒ EYEPOINT TECHNOLOGY®
- ☐ AS-WORN TECHNOLOGY™

NOTE:

- Prism to 10.00 Diopters

CR-39	TREATMENT						COATING
	Clear	Transitions Signature7™	XTRActive	Vantage	Polarized	Drivewear	
POWER RANGE	-14.00 to +8.00	-14.00 to +8.00	-14.00 to +8.00	-14.00 to +8.00	-14.00 to +8.00	-14.00 to +8.00	
ADD POWER							
CYLINDER	to -6.00	to -6.00	to -6.00	to -6.00	to -6.00	to -6.00	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15		
INVISIBLE MARKINGS	15 FF SV	TR FF SV	X5 FF SV	V5 FF SV	PZ FF SV	DW FF SV	

DLC™ TRIVEX™	TREATMENT					COATING
	Clear	Transitions Signature7™	XTRActive	Vantage	Polarized	
POWER RANGE	-14.00 to +7.50	-14.00 to +7.50	-14.00 to +7.50	-14.00 to +7.50	-14.00 to +7.50	
ADD POWER						
CYLINDER	to -6.50	to -6.50	to -6.50	to -6.50	to -6.50	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS	TX FF SV	TTX FF SV	XX FF SV	VX FF SV	PTX FF SV	

SuperLite™ 1.60	TREATMENT				COATING
	Clear	Transitions Signature7™	XTRActive	Polarized	
POWER RANGE	-14.00 to +8.00	-14.00 to +8.00	-14.00 to +8.00	-14.00 to +8.00	
ADD POWER					
CYLINDER	to -7.00	to -7.00	to -7.00	to -7.00	
OTHER		Gray/Brown	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS	16 FF SV	T6 FF SV	X6 FF SV	PZ6 FF SV	

1.56	TREATMENT		COATING
	Clear	Polarized	
POWER RANGE	-14.00 to +8.50	-14.00 to +8.50	
ADD POWER			
CYLINDER	to -6.75	to -6.75	
OTHER			
INVISIBLE MARKINGS	56 FF SV	56P FF SV	

SuperLite™ 1.67	TREATMENT				COATING
	Clear	Transitions Signature7™	XTRActive	Polarized	
POWER RANGE	-16.00 to +16.00	-16.00 to +16.00	-16.00 to +16.00	-16.00 to +16.00	
ADD POWER					
CYLINDER	to -8.00	to -8.00	to -8.00	to -8.00	
OTHER		Gray/Brown	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS	67 FF SV	T67 FF SV	X67 FF SV	PZ67 FF SV	

POLYPLUS™/ POLYCARBONATE	TREATMENT						COATING
	Clear	Transitions Signature7™	XTRActive	Vantage	Polarized	Drivewear	
POWER RANGE	-14.00 to +8.75	-14.00 to +8.75	-14.00 to +8.75	-14.00 to +8.75	-14.00 to +8.75	-14.00 to 8.75	
ADD POWER							
CYLINDER	to -7.00	to -7.00	to -7.00	to -7.00	to -7.00	to -7.00	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15		
INVISIBLE MARKINGS	PC FF SV	PT FF SV	PX FF SV	PV FF SV	PPZ FF SV	PD FF SV	

SuperLite™ 1.74	TREATMENT		COATING
	Clear	Transitions Signature7™	
POWER RANGE	-20.00 to +15.50	-16.00 to +12.50	
ADD POWER			
CYLINDER	to -9.00	to -9.00	
OTHER		Gray/Brown	
INVISIBLE MARKINGS	74 FF SV	T74 FF SV	

Crizal AR Coating
 Hoya AR Coating
 Zeiss AR Coating
 iCoat AR Coating

All Shamir products are available with your lab's in-house AR Coatings. Please check with your lab to determine which specific manufacturer's AR coating is available with the Shamir lens and material that you are ordering, as a manufacturer may offer more than one AR coating.




EFFECTIVE DATE: SEP. 15, 2014








- Denotes that material is available in a Shamir Freeform Blank™










EVERYDAY> SINGLE VISION

—SHAMIR
Single Vision™

CR-39	TREATMENT				
	Clear	Transitions	XTRActive	Vantage	Polarized
POWER RANGE	-10.00 to +6.00	-11.00 to +6.00	--	--	-8.00 to +6.00
ADD POWER					
CYLINDER	to -6.00	to -6.00	--	--	to -4.00
BASE CURVE	2, 3.50, 4.25, 6.25, 7.25	.5, 2.25, 4.25, 5.25, 6.25, 8.25	--	--	2.25, 3.25, 4.25, 5.25, 6.25, 7.25, 8.25
	Wrap: 4, 6, 7, 8				
OTHER/COATING	 Hard Coated/Uncoated	 Gray	--	--	 Gray/Brown
INVISIBLE MARKINGS					






SuperLite™ 1.60	TREATMENT				
	 Clear	 Transitions	 XTRActive	Vantage	Polarized
POWER RANGE	-15.00 to +9.00	-15.00 to +9.00	-15.00 to +9.00	--	-9.50 to +7.00
ADD POWER					
CYLINDER	to -7.00	to -7.00	to -7.00		to -4.00
BASE CURVE	.5, 2, 3, 4, 5, 6, 7, 8, 9	.5, 2, 3, 4, 5, 6, 7, 8, 9	.5, 2, 3, 4, 5, 6, 7, 8, 9	--	1.50, 2, 3, 4, 5, 6, 7, 8, 9
	Wrap: 6, 7, 8				
OTHER /COATING	 Hard Coated/Uncoated	 Gray/Brown	 Gray	--	 Gray/Brown
INVISIBLE MARKINGS					









SuperLite™ 1.67	TREATMENT				
	 Clear	 Transitions	 XTRActive	Vantage	Polarized
POWER RANGE	-16.00 to +16.25	-16.00 to +16.25	-16.00 to +16.25	--	-12.00 to +7.00
ADD POWER					
CYLINDER	to -8.00	to -8.00	to -8.00	--	to -4.00
BASE CURVE	.5, 2, 3, 4, 5, 6, 7, 8, 9, 11, 12.5	.5, 2, 3, 4, 5, 6, 7, 8, 9, 11, 12.5	.5, 2, 3, 4, 5, 6, 7, 8, 9, 11, 12.5	--	1.25, 2.75, 3.75, 4.75, 5.25, 6.75, 8
	Wrap: 6, 7, 8		Wrap: 6, 7, 8		
OTHER /COATING	 Hard Coated/Uncoated	 Gray/Brown	 Gray	--	 Gray/Brown
INVISIBLE MARKINGS					





PRODUCT FEATURES:

- ☐ FREEFORM PAL
- ☐ FREEFORM SV
- ☐ CONVENTIONAL PAL
- ☒ CONVENTIONAL SV
- ☐ EYEPOINT TECHNOLOGY®
- ☐ AS-WORN TECHNOLOGY®

SuperLite™ 1.74	TREATMENT				
	 Clear				
POWER RANGE	-18.00 to +8.50				
ADD POWER					
CYLINDER	to -9.00				
NOMINAL BASE CURVE	.5, 1, 2, 3, 4, 5, 6, 7, 8, 9				
	Extra Thick: .5, 1, 2				
OTHER/COATING	 Hard Coated/Uncoated				
INVISIBLE MARKINGS					


DLC™ TRIVEX™	TREATMENT			
	 Clear	 Transitions	XTRActive	Polarized
POWER RANGE	-14.00 to +7.50	-14.00 to +7.50	--	-7.00 to +7.50
ADD POWER				
CYLINDER	to -6.50	to -6.50	--	to -4.00
NOMINAL BASE CURVE	.5, 1, 2, 3, 4, 5, 6, 7, 8, 9 DLC™ Trivex™	.5, 1, 2, 3, 4, 5, 6, 7, 8, 9	--	2.25, 4.25, 6.25, 8.25
	Wrap: 6.25, 7.25, 8.25			
OTHER/COATING	 Hard Coated/Uncoated	 Gray/Brown	--	 Gray/Brown
INVISIBLE MARKINGS				

POLYPLUS™/ POLYCARBONATE	TREATMENT				
	 Clear	 Transitions	 XTRActive	Vantage	Polarized
POWER RANGE	-14.75 to +9.00	-14.75 to +9.00	-14.75 to +9.00	--	-10.00 to +6.00
ADD POWER					
CYLINDER	to -7.00	to -7.00	to -7.00	--	to -4.00
BASE CURVE	.5, 2, 3, 4, 5.25, 6.25, 7.5, 8.5, 9.5	.5, 2, 3, 4, 5.25, 6.25, 7.5, 8.5, 9.5	.5, 2, 3, 4, 5.25, 6.25, 7.5, 8.5, 9.5	--	1.50, 2.25, 3.25, 4.25, 5.50, 6.50, 7.50, 8.50
OTHER/COATING	Hard Coated:  Non Primer: 	 Gray/Brown	 Gray	--	 Gray/Brown
INVISIBLE MARKINGS					


 Crizal AR Coating  Hoya AR Coating  Zeiss AR Coating  iCoat AR Coating

All Shamir products are available with your lab's in-house AR Coatings. Please check with your lab to determine which specific manufacturer's AR coating is available with the Shamir lens and material that you are ordering, as a manufacturer may offer more than one AR coating.

EFFECTIVE DATE: MAY 9, 2013

 - Denotes that material is available in a Shamir Freeform Blank™


SHAMIR

NXT®	TREATMENT	
	Photochromic	Polarized
POWER RANGE	-5.00 to +4.00	-5.00 to +4.00
ADD POWER		
CYLINDER	to -6.00	to -4.00
BASE CURVE	6.00, 8.00	6.00, 8.00
OTHER/COATING	Green/Copper	 Gray/Brown
INVISIBLE MARKINGS		

PRODUCT FEATURES:

- ☐ FREEFORM PAL
- ☐ FREEFORM SV
- ☐ CONVENTIONAL PAL
- ☒ CONVENTIONAL SV
- ☐ EYEPOINT TECHNOLOGY®
- ☐ AS-WORN TECHNOLOGY®

PERFORMANCE

Shamir Attitude III® – Fashion

PRESCRIPTOR ORDER CODE:

ATTITUDE III FASHION

PRODUCT FEATURES:

- ☒ FREEFORM PAL
- ☐ FREEFORM OCCUPATION
- ☐ FREEFORM SV
- ☐ FREEFORM BIFOCAL
- ☐ CONVENTIONAL PAL
- ☐ CONVENTIONAL OCCUPATION
- ☐ CONVENTIONAL SV

- ☒ EYEPOINT TECHNOLOGY III®
- ☒ AS-WORN QUADRO™
- ☒ INTELICORRIDOR™
- ☒ NATURAL POSTURE™
- ☐ EYEPOINT TECHNOLOGY®
- ☐ AS-WORN TECHNOLOGY™

- NOTE:**
- Recommended min fitting heights are 15mm & 18mm
 - Prism to 10.00 Diopters

CR-39	TREATMENT						COATING
	Clear	Transitions Signature 7™	XTRActive	Vantage	Polarized	Drivewear	
POWER RANGE	-6.50 to +5.25	-6.50 to +5.25	-6.50 to +5.25	-6.50 to +5.25	-6.50 to +5.25	-6.50 to +5.25	
ADD POWER	0.75 to +3.50	0.75 to +3.50	0.75 to +3.50	0.75 to +3.50	0.75 to +3.50	0.75 to +3.00	
CYLINDER	to -6.00	to -6.00	to -6.00	to -6.00	to -6.00	to -6.00	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15		
INVISIBLE MARKINGS - 15mm	15 250	TR 250	X5 250	TV 250	PZ 250	DW 250	
INVISIBLE MARKINGS - 18mm	15 250	TR 250	X5 250	TV 250	PZ 250	DW 250	

SuperLite™ 1.60	TREATMENT					COATING
	Clear	Transitions Signature 7™	XTRActive	Polarized		
POWER RANGE	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00		
ADD POWER	0.75 to +3.50	0.75 to +3.50	0.75 to +3.50	0.75 to +3.50		
CYLINDER	to -7.00	to -7.00	to -7.00	to -7.00		
OTHER		Gray/Brown	Gray	Gray/Brown/G-15		
INVISIBLE MARKINGS - 15mm	16 250	T6 250	X6 250	PZ6 250		
INVISIBLE MARKINGS - 18mm	16 250	T6 250	X6 250	PZ6 250		

SuperLite™ 1.67	TREATMENT				COATING
	Clear	Transitions Signature 7™	XTRActive	Polarized	
POWER RANGE	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	
ADD POWER	0.75 to +3.50	0.75 to +3.50	0.75 to +3.50	0.75 to +3.50	
CYLINDER	to -8.00	to -8.00	to -8.00	to -8.00	
OTHER		Gray/Brown	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS - 15mm	67 250	T67 250	X67 250	PZ67 250	
INVISIBLE MARKINGS - 18mm	67 250	T67 250	X67 250	PZ67 250	

DLC™/TRIVEX™	TREATMENT					COATING
	Clear	Transitions Signature 7™	XTRActive	Vantage	Polarized	
POWER RANGE	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.00 to +6.00	
ADD POWER	0.75 to +3.50	0.75 to +3.50	0.75 to +3.50	0.75 to +3.50	0.75 to +3.50	
CYLINDER	to -6.50	to -6.50	to -6.50	to -6.50	to -6.50	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS - 15mm	TX 250	TTX 250	XX 250	VX 250	PTX 250	
INVISIBLE MARKINGS - 18mm	TX 250	TTX 250	XX 250	VX 250	PTX 250	

1.56	TREATMENT		COATING
	Clear	Polarized	
POWER RANGE	-6.50 to +6.00	-6.50 to +6.00	
ADD POWER	0.75 to +3.50	0.75 to +3.50	
CYLINDER	to -6.75	to -6.75	
OTHER	BluTech Indoor	BluTech Outdoor	
INVISIBLE MARKINGS - 15mm	56 250	56P 250	
INVISIBLE MARKINGS - 18mm	56 250	56P 250	

POLYPLUS™/ POLYCARBONATE	TREATMENT						COATING
	Clear	Transitions Signature 7™	XTRActive	Vantage	Polarized	Drivewear	
POWER RANGE	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	
ADD POWER	0.75 to +3.50	0.75 to +3.50	0.75 to +3.50	0.75 to +3.50	0.75 to +3.50	0.75 to +3.50	
CYLINDER	to -7.00	to -7.00	to -7.00	to -7.00	to -7.00	to -7.00	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15		
INVISIBLE MARKINGS - 15mm	PC 250	PT 250	PX 250	PV 250	PPZ 250	PD 250	
INVISIBLE MARKINGS - 18mm	PC 250	PT 250	PX 250	PV 250	PPZ 250	PD 250	

NXT	TREATMENT		COATING
	Photochromic	Polarized	
POWER RANGE	-6.00 to +5.50	-6.00 to +5.50	
ADD POWER	0.75 to +3.50	0.75 to +3.50	
CYLINDER	to -6.50	to -6.50	
OTHER	Green/Copper	Gray/Brown	
INVISIBLE MARKINGS - 15mm	NXTP 250	NXTZ 250	
INVISIBLE MARKINGS - 18mm	NXTP 250	NXTZ 250	

Crizal AR Coating Hoya AR Coating Zeiss AR Coating iCoat AR Coating

All Shamir products are available with your lab's in-house AR Coatings. Please check with your lab to determine which specific manufacturer's AR coating is available with the Shamir lens and material that you are ordering, as a manufacturer may offer more than one AR coating.

EFFECTIVE DATE: SEP. 15, 2014

- Denotes that material is available in a Shamir Freeform Blank™




































- ☒ FREEFORM PAL
- ☐ FREEFORM OCCUPATION
- ☐ FREEFORM SV
- ☐ FREEFORM BIFOCAL
- ☐ CONVENTIONAL PAL
- ☐ CONVENTIONAL OCCUPATION
- ☐ CONVENTIONAL SV
















- ☒ EYEPOINT TECHNOLOGY III®
- ☒ AS-WORN QUADRO™
- ☒ INTELLICORRIDOR™
- ☒ NATURAL POSTURE™
- ☐ EYEPOINT TECHNOLOGY®
- ☐ AS-WORN TECHNOLOGY™

























NOTE:
















- Recommended min fitting height is 18mm
- Prism to 10.00 Diopters









CR-39	TREATMENT						COATING
	Clear	Transitions Signature7™	XTRActive	Vantage	Polarized	Drivewear	
POWER RANGE	-6.50 to +5.25	-6.50 to +5.25	-6.50 to +5.25	-6.50 to +5.25	-6.50 to +5.25	-6.50 to +5.25	  
ADD POWER	0.75 to +3.50	0.75 to +3.50	0.75 to +3.50	0.75 to +3.50	0.75 to +3.50	0.75 to +3.50	
CYLINDER	to -6.00	to -6.00	to -6.00	to -6.00	to -6.00	to -6.00	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15		
INVISIBLE MARKINGS	 250  15  250	 250  TR  250	 250  X5  250	 250  TV  250	 250  PZ  250	 250  DW  250	



















1.56	TREATMENT		COATING
	Clear	Polarized	
POWER RANGE	-6.50 to +6.00	-6.50 to +6.00	  
ADD POWER	0.75 to +3.50	0.75 to +3.50	
CYLINDER	to -6.75	to -6.75	
OTHER	BluTech Indoor	BluTech Outdoor	
INVISIBLE MARKINGS	 250  56  250	 250  56P  250	  

SuperLite™ 1.60	TREATMENT				COATING
	Clear	Transitions Signature7™	XTRActive	Polarized	
POWER RANGE	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	  
ADD POWER	0.75 to +3.50	0.75 to +3.50	0.75 to +3.50	0.75 to +3.50	
CYLINDER	to -7.00	to -7.00	to -7.00	to -7.00	
OTHER		Gray/Brown	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS	 250  16  250	 250  T6  250	 250  X6  250	 250  PZ6  250	

POLYPLUS™/ POLYCARBONATE	TREATMENT						COATING
	Clear	Transitions Signature7™	XTRActive	Vantage	Polarized	Drivewear	
POWER RANGE	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	  
ADD POWER	0.75 to +3.50	0.75 to +3.50	0.75 to +3.50	0.75 to +3.50	0.75 to +3.50	0.75 to +3.50	
CYLINDER	to -7.00	to -7.00	to -7.00	to -7.00	to -7.00	to -7.00	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15		
INVISIBLE MARKINGS	 250  PC  250	 250  PT  250	 250  PX  250	 250  PV  250	 250  PPZ  250	 250  PD  250	  

SuperLite™ 1.67	TREATMENT				COATING
	Clear	Transitions Signature7™	XTRActive	Polarized	
POWER RANGE	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	  
ADD POWER	0.75 to +3.50	0.75 to +3.50	0.75 to +3.50	0.75 to +3.50	
CYLINDER	to -8.00	to -8.00	to -8.00	to -8.00	
OTHER		Gray/Brown	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS	 250  67  250	 250  T67  250	 250  X67  250	 250  PZ67  250	

NXT	TREATMENT		COATING
	Photochromic	Polarized	
POWER RANGE	-6.00 to +5.50	-6.00 to +5.50	
ADD POWER	0.75 to +3.50	0.75 to +3.50	
CYLINDER	to -6.50	to -6.50	
OTHER	Green/Copper	Gray/Brown	
INVISIBLE MARKINGS	 250  NXTP  250	 250  NXTZ  250	

DLC™/TRIVEX™	TREATMENT					COATING
	Clear	Transitions Signature7™	XTRActive	Vantage	Polarized	
POWER RANGE	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	  
ADD POWER	0.75 to +3.50	0.75 to +3.50	0.75 to +3.50	0.75 to +3.50	0.75 to +3.50	
CYLINDER	to -6.50	to -6.50	to -6.50	to -6.50	to -6.50	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS	 250  TX  250	 250  TTX  250	 250  XX  250	 250  VX  250	 250  PTX  250	



PREScriPTOR ORDER CODE:
GOLF

- PRODUCT FEATURES:
- ☒ FREEFORM PAL
 - ☐ FREEFORM SV
 - ☐ CONVENTIONAL PAL
 - ☐ CONVENTIONAL SV
 - ☒ EYEPOINT TECHNOLOGY®
 - ☒ AS-WORN TECHNOLOGY®

NOTE:

- Recommended min fitting height is 18mm
- Prism to 10.00 Diopters

SuperLite™ 1.60	TREATMENT				COATING
	Clear	Transitions® Signature7™	XTRActive	Polarized	
POWER RANGE	-6.00 to +5.25	-6.00 to +5.25	-6.00 to +5.25	-6.00 to +5.25	
ADD POWER	+1.00 to +3.50	+1.00 to +3.50	+1.00 to +3.50	+1.00 to +3.50	
CYLINDER	to -7.00	to -7.00	to -7.00	to -7.00	
OTHER		Gray/Brown	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS	350 16	350 T6	350 X6	350 PZ6	

NXT	TREATMENT		COATING
	Photochromic	Polarized	
POWER RANGE	-5.00 to +4.00	-5.00 to +4.00	
ADD POWER	+1.00 to +3.50	+1.00 to +3.50	
CYLINDER	to -6.50	to -6.50	
OTHER			
INVISIBLE MARKINGS	350 NXTP	350 NXTZ	

DLC™ TRIVEX™	TREATMENT					COATING
	Clear	Transitions® Signature7™	XTRActive	Vantage	Polarized	
POWER RANGE	-5.00 to +4.00	-5.00 to +4.00	-5.00 to +4.00	-5.00 to +4.00	-5.00 to +4.00	
ADD POWER	+1.00 to +3.50	+1.00 to +3.50	+1.00 to +3.50	+1.00 to +3.50	+1.00 to +3.50	
CYLINDER	to -6.50	to -6.50	to -6.50	to -6.50	to -6.50	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS	350 TX	350 TTX	350 XX	350 VX	350 PTX	

1.56	TREATMENT		COATING
	Clear	Polarized	
POWER RANGE	-5.50 to +5.25	-5.50 to +5.25	
ADD POWER	+1.00 to +3.50	+1.00 to +3.50	
CYLINDER	to -6.75	to -6.75	
OTHER			
INVISIBLE MARKINGS	350 56	350 56P	

POLYPLUS™/ POLYCARBONATE	TREATMENT						COATING
	Clear	Transitions® Signature7™	XTRActive	Vantage	Polarized	Drivewear	
POWER RANGE	-6.00 to +5.25	-6.00 to +5.25	-6.00 to +5.25	-6.00 to +5.25	-6.00 to +5.25	-6.00 to +5.25	
ADD POWER	+1.00 to +3.50	+1.00 to +3.50	+1.00 to +3.50	+1.00 to +3.50	+1.00 to +3.50	+1.00 to +3.50	
CYLINDER	to -7.00	to -7.00	to -7.00	to -7.00	to -7.00	to -7.00	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15		
INVISIBLE MARKINGS	350 PC	350 PT	350 PX	350 PPT	350 PPZ	350 PD	

Crizal AR Coating Hoya AR Coating Zeiss AR Coating iCoat AR Coating

All Shamir products are available with your lab's in-house AR Coatings. Please check with your lab to determine which specific manufacturer's AR coating is available with the Shamir lens and material that you are ordering, as a manufacturer may offer more than one AR coating.

EFFECTIVE DATE: SEP. 15, 2014

- Denotes that material is available in a Shamir Freeform Blank™



Autograph II

ATTITUDE™

PRESCRIPTION ORDER CODE:

AUTOGRAPH II ATTITUDE 15
AUTOGRAPH II ATTITUDE 18

PRODUCT FEATURES:

- ☒ FREEFORM PAL
- ☐ FREEFORM SV
- ☐ CONVENTIONAL PAL
- ☐ CONVENTIONAL SV
- ☒ EYEPOINT TECHNOLOGY®
- ☒ AS-WORN TECHNOLOGY™

NOTE:

- Recommended min fitting heights are 15mm and 18mm
- Prism to 10.00 Diopters

CR-39	TREATMENT						COATING
	Clear	Transitions Signature7™	XTRActive	Vantage	Polarized	Drivewear	
POWER RANGE	-5.00 to +4.00	-5.00 to +4.00	-5.00 to +4.00	-5.00 to +4.00	-5.00 to +4.00	-5.00 to +4.00	
ADD POWER	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	
CYLINDER	to -6.00	to -6.00	to -6.00	to -6.00	to -6.00	to -6.00	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15		
INVISIBLE MARKINGS - 15mm	15 250 FF-AS	TR 250 FF-AS	X5 250 FF-AS	V5 250 FF-AS	PZ 250 FF-AS	DW 250 FF-AS	
INVISIBLE MARKINGS - 18mm	15 250 FFA	TR 250 FFA	X5 250 FFA	V5 250 FFA	PZ 250 FFA	DW 250 FFA	

DLC™ TRIVEX™	TREATMENT					COATING
	Clear	Transitions Signature7™	XTRActive	Vantage	Polarized	
POWER RANGE	-5.00 to +4.00	-5.00 to +4.00	-5.00 to +4.00	-5.00 to +4.00	-5.00 to +4.00	
ADD POWER	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	
CYLINDER	to -6.50	to -6.50	to -6.50	to -6.50	to -6.50	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS - 15mm	TX 250 FF-AS	TTX 250 FF-AS	XX 250 FF-AS	VX 250 FF-AS	PTX 250 FF-AS	
INVISIBLE MARKINGS - 18mm	TX 250 FFA	TTX 250 FFA	XX 250 FFA	VX 250 FFA	PTX 250 FFA	

SuperLite™ 1.60	TREATMENT				COATING
	Clear	Transitions Signature7™	XTRActive	Polarized	
POWER RANGE	-6.50 to +6.00	-6.50 to +6.00	-6.00 to +5.25	-6.00 to +5.25	
ADD POWER	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	
CYLINDER	to -7.00	to -7.00	to -7.00	to -7.00	
OTHER		Gray/Brown	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS - 15mm	16 250 FF-AS	T6 250 FF-AS	X6 250 FF-AS	PZ6 250 FF-AS	
INVISIBLE MARKINGS - 18mm	16 250 FFA	T6 250 FFA	X6 250 FFA	PZ6 250 FFA	

POLYPLUS™/ POLYCARBONATE	TREATMENT						COATING
	Clear	Transitions Signature7™	XTRActive	Vantage	Polarized	Drivewear	
POWER RANGE	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	
ADD POWER	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	
CYLINDER	to -7.00	to -7.00	to -7.00	to -7.00	to -7.00	to -7.00	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15		
INVISIBLE MARKINGS - 15mm	PC 250 FF-AS	PT 250 FF-AS	PX 250 FF-AS	PV 250 FF-AS	PPZ 250 FF-AS	PD 250 FF-AS	
INVISIBLE MARKINGS - 18mm	PC 250 FFA	PT 250 FFA	PX 250 FFA	PV 250 FFA	PPZ 250 FFA	PD 250 FFA	

SuperLite™ 1.67	TREATMENT				COATING
	Clear	Transitions Signature7™	XTRActive	Polarized	
POWER RANGE	-6.00 to +5.00	-6.00 to +5.00	-6.00 to +5.00	-6.00 to +5.00	
ADD POWER	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	
CYLINDER	to -8.00	to -8.00	to -8.00	to -8.00	
OTHER		Gray/Brown	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS - 15mm	67 250 FF-AS	T67 250 FF-AS	X67 250 FF-AS	PZ67 250 FF-AS	
INVISIBLE MARKINGS - 18mm	67 250 FFA	T67 250 FFA	X67 250 FFA	PZ67 250 FFA	

NXT	TREATMENT		COATING
	Photochromic	Polarized	
POWER RANGE	-5.00 to +4.00	-5.00 to +4.00	
ADD POWER	+0.75 to +3.50	+0.75 to +3.50	
CYLINDER	to -6.50	to -6.50	
OTHER	Green/Copper	Gray/Brown	
INVISIBLE MARKINGS - 15mm	NXTP 250 FF-AS	NXTZ 250 FF-AS	
INVISIBLE MARKINGS - 18mm	NXTP 250 FFA	NXTZ 250 FFA	

PERFORMANCE> SPORT

SHAMIR

Autograph II

ATTITUDE™

PREScriPTOR ORDER CODE:

AUTOGRAPH II ATTITUDE 15

AUTOGRAPH II ATTITUDE 18

- PRODUCT FEATURES:
- ☒ FREEFORM PAL

☐ FREEFORM SV

☐ CONVENTIONAL PAL








☐ CONVENTIONAL SV
- ☒ EYEPOINT TECHNOLOGY®

☒ AS-WORN TECHNOLOGY™

NOTE:

• Recommended min fitting heights are 15mm and 18mm

• Prism to 10.00 Diopters

1.56	TREATMENT		COATING
	Clear	Polarized	
POWER RANGE	-5.50 to +6.00	-5.50 to +6.00	
ADD POWER	+0.75 to +3.50	+0.75 to +3.50	
CYLINDER	to -6.75	to -6.75	
OTHER			 
INVISIBLE MARKINGS - 15mm	 56 250 FF-AS	 56P 250 FF-AS	
INVISIBLE MARKINGS - 18mm	 56 250 FFA	 56P 250 FFA	

CR-39	TREATMENT				
	Clear	Transitions®	XTRActive	Vantage	Polarized
POWER RANGE	--	--	--	--	--
ADD POWER	--	--	--	--	--
BASE CURVE	--	--	--	--	--
OTHER/COATING	--	--	--	--	--
INVISIBLE MARKINGS	--	--	--	--	--

SuperLite™ 1.60	TREATMENT				
	Clear	Transitions®	XTRActive	Vantage	Polarized
POWER RANGE	--	--	--	--	--
ADD POWER	--	--	--	--	--
BASE CURVE	--	--	--	--	--
OTHER/COATING	--	--	--	--	--
INVISIBLE MARKINGS	--	--	--	--	--

SuperLite™ 1.67	TREATMENT				
	Clear	Transitions®	XTRActive	Vantage	Polarized
POWER RANGE	--	--	--	--	--
ADD POWER	--	--	--	--	--
BASE CURVE	--	--	--	--	--
OTHER/COATING	--	--	--	--	--
INVISIBLE MARKINGS	--	--	--	--	--

PRODUCT FEATURES:





- ☐ FREEFORM PAL
- ☐ FREEFORM SV
- ☒ CONVENTIONAL PAL
- ☐ CONVENTIONAL SV
- ☒ EYEPOINT TECHNOLOGY®
- ☐ AS-WORN TECHNOLOGY®

NOTE:

- Recommended min fitting height is 16mm
- Cylinder to -6.00
- Prism to 4.00 Diopters

SuperLite™ 1.74	TREATMENT				
	Clear				
POWER RANGE	--				
ADD POWER	--				
BASE CURVE	--				
OTHER/COATING	--				
INVISIBLE MARKINGS	--				

DLC™ TRIVEX™	TREATMENT			
	Clear	Transitions®	XTRActive	Polarized
POWER RANGE	--	--	--	--
ADD POWER	--	--	--	--
BASE CURVE	--	--	--	--
OTHER/COATING	--	--	--	--
INVISIBLE MARKINGS	--	--	--	--

POLYPLUS™ / POLYCARBONATE	TREATMENT				
	Clear	Transitions®	XTRActive	Vantage	Polarized
POWER RANGE	-4.00 to +2.00	--	--	--	-4.00 to +2.00
ADD POWER	+1.00 to +3.00	--	--	--	+1.00 to +3.00
BASE CURVE	8	--	--	--	8
OTHER/COATING	 Piccolo®	--	--	--	 Piccolo® Gray/Brown
INVISIBLE MARKINGS	 250	--	--	--	 250

PERFORMANCE> SPORT SV

Shamir Attitude III® – SV

PRESCRIPTOR ORDER CODE:

ATTITUDE III SV

PRODUCT FEATURES:

- ☐ FREEFORM PAL
- ☐ FREEFORM OCCUPATION
- ☒ FREEFORM SV
- ☐ FREEFORM BIFOCAL
- ☐ CONVENTIONAL PAL
- ☐ CONVENTIONAL OCCUPATION
- ☐ CONVENTIONAL SV
- ☒ EYEPOINT TECHNOLOGY III®
- ☒ AS-WORN QUADRO™
- ☐ INTELLICORRIDOR™
- ☐ NATURAL POSTURE™
- ☐ EYEPOINT TECHNOLOGY®
- ☐ AS-WORN TECHNOLOGY™

NOTE:

- Prism to 10.00 Diopters

CR-39	TREATMENT						COATING
	Clear	Transitions Signature7™	XTRActive	Vantage	Polarized	Drivewear	
POWER RANGE	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	
CYLINDER	to -6.00	to -6.00	to -6.00	to -6.00	to -6.00	to -6.00	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15		
INVISIBLE MARKINGS	15 250 AT ₃	TR 250 AT ₃	X5 250 AT ₃	TV 250 AT ₃	PZ 250 AT ₃	DW 250 AT ₃	

1.56	TREATMENT		COATING
	Clear	Polarized	
POWER RANGE	-6.50 to +6.00	-6.50 to +6.00	
CYLINDER	to -6.75	to -6.75	
OTHER	BluTech Indoor	BluTech Outdoor	
INVISIBLE MARKINGS	56 250 AT ₃	56P 250 AT ₃	

SuperLite™ 1.60	TREATMENT				COATING
	Clear	Transitions Signature7™	XTRActive	Polarized	
POWER RANGE	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	
CYLINDER	to -7.00	to -7.00	to -7.00	to -7.00	
OTHER		Gray/Brown	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS	16 250 AT ₃	T6 250 AT ₃	X6 250 AT ₃	PZ6 250 AT ₃	

POLYPLUS™/ POLYCARBONATE	TREATMENT						COATING
	Clear	Transitions Signature7™	XTRActive	Vantage	Polarized	Drivewear	
POWER RANGE	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	
CYLINDER	to -7.00	to -7.00	to -7.00	to -7.00	to -7.00	to -7.00	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15		
INVISIBLE MARKINGS	PC 250 AT ₃	PT 250 AT ₃	PX 250 AT ₃	PV 250 AT ₃	PPZ 250 AT ₃	PD 250 AT ₃	

SuperLite™ 1.67	TREATMENT				COATING
	Clear	Transitions Signature7™	XTRActive	Polarized	
POWER RANGE	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	
CYLINDER	to -8.00	to -8.00	to -8.00	to -8.00	
OTHER		Gray/Brown	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS	250 AT ₃	T67 250 AT ₃	X67 250 AT ₃	PZ67 250 AT ₃	

NXT	TREATMENT		COATING
	Photochromic	Polarized	
POWER RANGE	-6.00 to +6.00	-6.00 to +6.00	
CYLINDER	to -6.50	to -6.50	
OTHER	Green/Copper	Gray/Brown	
INVISIBLE MARKINGS	NXT 250 AT ₃	NXTZ 250 AT ₃	

DLC™/TRIVEX™	TREATMENT					COATING
	Clear	Transitions Signature7™	XTRActive	Vantage	Polarized	
POWER RANGE	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	
CYLINDER	to -6.50	to -6.50	to -6.50	to -6.50	to -6.50	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS	TX 250 AT ₃	TTX 250 AT ₃	XX 250 AT ₃	VX 250 AT ₃	PTX 250 AT ₃	

Crizal AR Coating Hoya AR Coating Zeiss AR Coating iCoat AR Coating

All Shamir products are available with your lab's in-house AR Coatings. Please check with your lab to determine which specific manufacturer's AR coating is available with the Shamir lens and material that you are ordering, as a manufacturer may offer more than one AR coating.

EFFECTIVE DATE: SEP. 15, 2014

- Denotes that material is available in a Shamir Freeform Blank™



PERFORMANCE> SPORT SV



PREScriPTOR ORDER CODE:
GOLF SV

- PRODUCT FEATURES:
- ☐ FREEFORM PAL
 - ☒ FREEFORM SV
 - ☐ CONVENTIONAL PAL
 - ☐ CONVENTIONAL SV
 - ☒ EYEPOINT TECHNOLOGY®
 - ☒ AS-WORN TECHNOLOGY™

NOTE:
• Prism to 10.00 Diopters

SuperLite™ 1.60	TREATMENT				COATING
	Clear	Transitions Signature7™	XTRActive	Polarized	
	POWER RANGE	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	
	CYLINDER	to -7.00	to -7.00	to -7.00	
	OTHER		Gray/Brown	Gray	
INVISIBLE MARKINGS	16 SV	T6 SV	X6 SV	PZ6 SV	

NXT	TREATMENT		COATING
	Photochromic	Polarized	
	POWER RANGE	-6.00 to +6.00	
	CYLINDER	to -6.50	
	OTHER		
INVISIBLE MARKINGS	NXTP SV	NXTZ SV	

DLC™ TRIVEX™	TREATMENT				COATING
	Clear	Transitions Signature7™	XTRActive	Polarized	
	POWER RANGE	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	
	CYLINDER	to -6.50	to -6.50	to -6.50	
	OTHER		Gray/Brown	Gray	
INVISIBLE MARKINGS	TX SV	TTX SV	XX SV	PTX SV	

1.56	TREATMENT		COATING
	Clear	Polarized	
	POWER RANGE	-6.00 to +6.00	
	CYLINDER	to -6.75	
	OTHER		
INVISIBLE MARKINGS	56 SV	56P SV	

POLYPLUS™/ POLYCARBONATE	TREATMENT						COATING
	Clear	Transitions Signature7™	XTRActive	Vantage	Polarized	Drivewear	
	POWER RANGE	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	
	CYLINDER	to -7.00	to -7.00	to -7.00	to -7.00	to -7.00	
	OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS	PC SV	PT SV	PX SV	PPT SV	PPZ SV	PD SV	

Autograph II

SINGLE VISION ATTITUDE™

PRESCRIBTOR ORDER CODE:



AUTOGRAPH II SV ATTITUDE





PRODUCT FEATURES:





- ☐ FREEFORM PAL
- ☒ FREEFORM SV
- ☐ CONVENTIONAL PAL
- ☐ CONVENTIONAL SV
- ☒ EYEPOINT TECHNOLOGY®
- ☒ AS-WORN TECHNOLOGY™





NOTE:





• Prism to 10.00 Diopters


CR-39	TREATMENT						COATING
	Clear	Transitions Signature7™	XTRActive	Vantage	Polarized	Drivewear	
POWER RANGE	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	   
ADD POWER							
CYLINDER	to -6.00	to -6.00	to -6.00	to -6.00	to -6.00	to -6.00	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15		
INVISIBLE MARKINGS	FF SA-15	FF SA-TR	FF SA-X5	FF SA-V5	FF SA-PZ	FF SA-DW	

DLC™ TRIVEX™	TREATMENT					COATING
	Clear	Transitions Signature7™	XTRActive	Vantage	Polarized	
POWER RANGE	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	   
ADD POWER						
CYLINDER	to -6.50	to -6.50	to -6.50	to -6.50	to -6.50	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS	FF SA-TX	FF SA-TTX	FF SA-XX	FF SA-VX	FF SA-PTX	

SuperLite™ 1.60	TREATMENT				COATING
	Clear	Transitions Signature7™	XTRActive	Polarized	
POWER RANGE	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	   
ADD POWER					
CYLINDER	to -7.00	to -7.00	to -7.00	to -7.00	
OTHER		Gray/Brown	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS	FF SA-16	FF SA-T6	FF SA-X6	FF SA-PZ6	

POLYPLUS™ / POLYCARBONATE	TREATMENT						COATING
	Clear	Transitions Signature7™	XTRActive	Vantage	Polarized	Drivewear	
POWER RANGE	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	-6.50 to +6.00	   
ADD POWER							
CYLINDER	to -7.00	to -7.00	to -7.00	to -7.00	to -7.00	to -7.00	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15		
INVISIBLE MARKINGS	FF SA-PC	FF SA-PT	FF SA-PX	FF SA-PV	FF SA-PPZ	FF SA-PD	

SuperLite™ 1.67	TREATMENT				COATING
	Clear	Transitions Signature7™	XTRActive	Polarized	
POWER RANGE	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	   
ADD POWER					
CYLINDER	to -8.00	to -8.00	to -8.00	to -8.00	
OTHER		Gray/Brown	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS	FF SA-67	FF SA-T67	FF SA-X67	FF SA-PZ67	

NXT	TREATMENT		COATING
	Photochromic	Polarized	
POWER RANGE	-6.00 to +6.00	-6.00 to +6.00	
ADD POWER			
CYLINDER	to -6.50	to -6.50	
OTHER	Green/Copper	Gray/Brown	
INVISIBLE MARKINGS	250 SA-NXTP	250 SA-NXTZ	

Autograph II

SINGLE VISION ATTITUDE™

PREScriPTOR ORDER CODE:

AUTOGRAPH II SV ATTITUDE

PRODUCT FEATURES:

- ☐ FREEFORM PAL
- ☒ FREEFORM SV
- ☐ CONVENTIONAL PAL
- ☐ CONVENTIONAL SV
- ☒ EYEPOINT TECHNOLOGY®
- ☒ AS-WORN TECHNOLOGY™

NOTE:

- Prism to 10.00 Diopters

1.56	TREATMENT		COATING
	Clear	Polarized	
POWER RANGE	-6.00 to +6.00	-6.00 to +6.00	<div><div></div><div></div><div></div><div></div></div>
ADD POWER			
CYLINDER	to -6.75	to -6.75	
OTHER			
INVISIBLE MARKINGS	<div><div></div><div></div></div> 250 SA-56	<div><div></div><div></div></div> 250 SA-56P	

OCCUPATION

Shamir Computer™

PRESCRITOR ORDER CODE:

COMPUTER

PRODUCT FEATURES:

- ☐ FREEFORM PAL
☒ FREEFORM OCCUPATION
☐ FREEFORM SV
☐ CONVENTIONAL PAL
☐ CONVENTIONAL OCCUPATION
☐ CONVENTIONAL SV
- ☐ EYEPOINT TECHNOLOGY III®
☐ AS-WORN QUADRO™
☒ INTELICORRIDOR™
☒ NATURAL POSTURE™
- ☐ EYEPOINT TECHNOLOGY®
☐ AS-WORN TECHNOLOGY™

NOTE:

- Recommended min fitting height is 14mm lower and 10mm upper
- Prism to 10.00 Diopters

CR-39	TREATMENT				COATING
	Clear	Transitions Signature7™	XTRActive	Polarized	
TOTAL NEAR POWER RANGE	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	
DYNAMIC REDUCTION	-0.75 to -2.00	-0.75 to -2.00	-0.75 to -2.00	-0.75 to -2.00	
CYLINDER	to -6.00	to -6.00	to -6.00	to -6.00	
OTHER		Gray/Brown	Gray	Gray/Brown/G-15	
INVISIBLE MARKING	15 250 FFC	TR 250 FFC	X5 250 FFC	PP2 250 FFC	

DLC™/TRIVEX™	TREATMENT			COATING
	Clear	Transitions Signature7™	XTRActive	
TOTAL NEAR POWER RANGE	-6.00 to +8.00	-6.00 to +8.00	-6.00 to +8.00	
DYNAMIC REDUCTION	-0.75 to -2.00	-0.75 to -2.00	-0.75 to -2.00	
CYLINDER	to -6.50	to -6.50	to -6.50	
OTHER		Gray/Brown	Gray	
INVISIBLE MARKING	TX 250 FFC	TTX 250 FFC	XX 250 FFC	

SuperLite™ 1.60	TREATMENT			COATING
	Clear	Transitions Signature7™	XTRActive	
TOTAL NEAR POWER RANGE	-6.00 to +8.00	-6.00 to +8.00	-6.00 to +8.00	
DYNAMIC REDUCTION	-0.75 to -2.00	-0.75 to -2.00	-0.75 to -2.00	
CYLINDER	to -7.00	to -7.00	to -7.00	
OTHER		Gray/Brown	Gray	
INVISIBLE MARKING	16 250 FFC	T6 250 FFC	X6 250 FFC	

1.56	TREATMENT		COATING
	Clear	Polarized	
TOTAL NEAR POWER RANGE	-6.00 to +8.00	-6.00 to +8.00	
DYNAMIC REDUCTION	-0.75 to -2.00	-0.75 to -2.00	
CYLINDER	to -6.75	to -6.75	
OTHER	BluTech Indoor	BluTech Indoor	
INVISIBLE MARKING	56 250 FFC	56P 250 FFC	

POLYPLUS™/ POLYCARBONATE	TREATMENT			COATING
	Clear	Transitions Signature7™	XTRActive	
TOTAL NEAR POWER RANGE	-6.00 to +8.00	-6.00 to +8.00	-6.00 to +8.00	
DYNAMIC REDUCTION	-0.75 to -2.00	-0.75 to -2.00	-0.75 to -2.00	
CYLINDER	to -7.00	to -7.00	to -7.00	
OTHER		Gray/Brown	Gray	
INVISIBLE MARKING	PC 250 FFC	PT 250 FFC	PX 250 FFC	

Crizal AR Coating
 Hoya AR Coating
 Zeiss AR Coating
 iCoat AR Coating

All Shamir products are available with your lab's in-house AR Coatings. Please check with your lab to determine which specific manufacturer's AR coating is available with the Shamir lens and material that you are ordering, as a manufacturer may offer more than one AR coating.

- Denotes that material is available in a Shamir Freeform Blank™

EFFECTIVE DATE: SEP. 15, 2014

OCCUPATION

Shamir WorkSpace™

PRESCRIBTOR ORDER CODE:

WORKSPACE

PRODUCT FEATURES:

- ☐ FREEFORM PAL
☒ FREEFORM OCCUPATION
☐ FREEFORM SV
☐ CONVENTIONAL PAL
☐ CONVENTIONAL OCCUPATION
☐ CONVENTIONAL SV
- ☐ EYEPOINT TECHNOLOGY III®
☐ AS-WORN QUADRO™
☒ INTELICORRIDOR™
☒ NATURAL POSTURE™
- ☐ EYEPOINT TECHNOLOGY®
☐ AS-WORN TECHNOLOGY™

NOTE:

- Recommended min fitting height is 14mm lower and 10mm upper
- Prism to 10.00 Diopters

CR-39	TREATMENT				COATING
	Clear	Transitions Signature7™	XTRActive	Polarized	
TOTAL NEAR POWER RANGE	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	-6.00 to +6.00	
DYNAMIC REDUCTION	-0.75 to -2.25	-0.75 to -2.50	-0.75 to -2.50	-0.75 to -2.50	
CYLINDER	to -6.00	to -6.00	to -6.00	to -6.00	
OTHER		Gray/Brown	Gray	Gray/Brown/G-15	
INVISIBLE MARKING	15 250 FFW	TR 250 FFW	X5 250 FFW	PPZ 250 FFW	

DLC™/TRIVEX™	TREATMENT			COATING
	Clear	Transitions Signature7™	XTRActive	
TOTAL NEAR POWER RANGE	-6.00 to +8.00	-6.00 to +8.00	-6.00 to +8.00	
DYNAMIC REDUCTION	-0.75 to -2.50	-0.75 to -2.50	-0.75 to -2.50	
CYLINDER	to -6.50	to -6.50	to -6.50	
OTHER		Gray/Brown	Gray	
INVISIBLE MARKING	TX 250 FFW	TTX 250 FFW	XX 250 FFW	

SuperLite™ 1.60	TREATMENT			COATING
	Clear	Transitions Signature7™	XTRActive	
TOTAL NEAR POWER RANGE	-6.00 to +8.00	-6.00 to +8.00	-6.00 to +8.00	
DYNAMIC REDUCTION	-0.75 to -2.50	-0.75 to -2.50	-0.75 to -2.50	
CYLINDER	to -7.00	to -7.00	to -7.00	
OTHER		Gray/Brown	Gray	
INVISIBLE MARKING	16 250 FFW	T6 250 FFW	X6 250 FFW	

1.56	TREATMENT		COATING
	Clear	Polarized	
TOTAL NEAR POWER RANGE	-6.00 to +8.00	-6.00 to +8.00	
DYNAMIC REDUCTION	-0.75 to -2.50	-0.75 to -2.50	
CYLINDER	to -6.75	to -6.75	
OTHER	BlueTech Indoor	BlueTech Outdoor	
INVISIBLE MARKING	56 250 FFW	56P 250 FFW	

POLYPLUS™/ POLYCARBONATE	TREATMENT			COATING
	Clear	Transitions Signature7™	XTRActive	
TOTAL NEAR POWER RANGE	-6.00 to +8.00	-6.00 to +8.00	-6.00 to +8.00	
DYNAMIC REDUCTION	-0.75 to -2.50	-0.75 to -2.50	-0.75 to -2.50	
CYLINDER	to -7.00	to -7.00	to -7.00	
OTHER		Gray/Brown	Gray	
INVISIBLE MARKING	PC 250 FFW	PT 250 FFW	PX 250 FFW	

Crizal AR Coating
 Hoya AR Coating
 Zeiss AR Coating
 iCoat AR Coating

All Shamir products are available with your lab's in-house AR Coatings. Please check with your lab to determine which specific manufacturer's AR coating is available with the Shamir lens and material that you are ordering, as a manufacturer may offer more than one AR coating.

- Denotes that material is available in a Shamir Freeform Blank™

EFFECTIVE DATE: SEP. 15, 2014

OCCUPATION

SHAMIR

Autograph II

OFFICE™

PRESCRIPTOR ORDER CODE:








AUTOGRAPH II OFFICE








PRODUCT FEATURES:








- ☐ FREEFORM PAL
- ☒ FREEFORM OCCUPATION
- ☐ FREEFORM SV
- ☐ CONVENTIONAL PAL
- ☐ CONVENTIONAL OCCUPATION
- ☐ CONVENTIONAL SV
- ☒ EYEPOINT TECHNOLOGY®
- ☐ AS-WORN TECHNOLOGY®





NOTE:

- Recommended min fitting height is 16mm
- Prism to 10.00 Diopters

CR-39	TREATMENT	COATING
	 Clear	
	TOTAL NEAR POWER RANGE	-9.50 to +5.25
	DYNAMIC REDUCTION	-0.75 to -2.25
	CYLINDER	to -6.00
	OTHER	
INVISIBLE MARKINGS	 -175  15 FFO	   


SuperLite™ 1.60	TREATMENT	COATING
	 Clear	
	TOTAL NEAR POWER RANGE	-9.50 to +7.25
	DYNAMIC REDUCTION	-0.75 to -2.25
	CYLINDER	to -7.00
	OTHER	
INVISIBLE MARKINGS	 -175  16 FFO	   

POLYPLUS™/ POLYCARBONATE	TREATMENT	COATING
	 Clear	
	TOTAL NEAR POWER RANGE	-9.50 to +7.25
	DYNAMIC REDUCTION	-0.75 to -2.25
	CYLINDER	to -7.00
	OTHER	
INVISIBLE MARKINGS	 -175  PC FFO	   

-  Crizal AR Coating
-  Hoya AR Coating
-  Zeiss AR Coating
-  iCoat AR Coating

All Shamir products are available with your lab's in-house AR Coatings. Please check with your lab to determine which specific manufacturer's AR coating is available with the Shamir lens and material that you are ordering, as a manufacturer may offer more than one AR coating.

EFFECTIVE DATE: SEP. 15, 2014

 - Denotes that material is available in a Shamir Freeform Blank™



OCCUPATION



PRESCRIPTOR ORDER CODE:

RELAX

PRODUCT FEATURES:

- ☐ FREEFORM PAL
- ☒ FREEFORM OCCUPATION
- ☒ FREEFORM SV
- ☐ CONVENTIONAL PAL
- ☐ CONVENTIONAL SV
- ☐ CONVENTIONAL OCCUPATION
- ☒ EYEPOINT TECHNOLOGY™
- ☐ AS-WORN TECHNOLOGY™

- NOTE:
- Recommended min fitting height is 16mm

• Fatigue relief solution with .65 Add in lower part of lens

• Prism to 10.00 Diopters

CR-39	TREATMENT						COATING
	Clear	Transitions Signature 7™	XTRActive	Vantage	Polarized	Drivewear	
POWER RANGE	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	
ADD POWER							
CYLINDER	to -6.00	to -6.00	to -6.00	to -6.00	to -6.00	to -6.00	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15		
INVISIBLE MARKINGS	15 RLX	TR RLX	X5 RLX	V5 RLX	PZ RLX	DW RLX	

DLC™ TRIVEX™	TREATMENT				COATING
	Clear	Transitions Signature 7™	XTRActive	Polarized	
POWER RANGE	-13.25 to +6.00	-13.25 to +6.00	-13.25 to +6.00	-13.25 to +6.00	
ADD POWER					
CYLINDER	to -6.50	to -6.50	to -6.50	to -6.50	
OTHER		Gray/Brown	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS	TX RLX	TTX RLX	XX RLX	PTX RLX	

SuperLite™ 1.60	TREATMENT				COATING
	Clear	Transitions Signature 7™	XTRActive	Polarized	
POWER RANGE	-14.00 to +8.00	-14.00 to +8.00	-14.00 to +8.00	-14.00 to +8.00	
ADD POWER					
CYLINDER	to -7.00	to -7.00	to -7.00	to -7.00	
OTHER		Gray/Brown	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS	16 RLX	T6 RLX	X6 RLX	PZ6 RLX	

1.56	TREATMENT		COATING
	Clear	Polarized	
POWER RANGE	-14.00 to +6.50	-14.00 to +6.50	
ADD POWER			
CYLINDER	to -6.75	to -6.75	
OTHER			
INVISIBLE MARKINGS	56 RLX	56P RLX	

SuperLite™ 1.67	TREATMENT				COATING
	Clear	Transitions Signature 7™	XTRActive	Polarized	
POWER RANGE	-16.75 to +8.25	-16.75 to +8.25	-16.75 to +8.25	-16.75 to +8.25	
ADD POWER					
CYLINDER	to -8.00	to -8.00	to -8.00	to -8.00	
OTHER		Gray/Brown	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS	67 RLX	T67 RLX	X67 RLX	PZ67 RLX	

POLYPLUS™/ POLYCARBONATE	TREATMENT						COATING
	Clear	Transitions Signature 7™	XTRActive	Vantage	Polarized	Drivewear	
POWER RANGE	-14.00 to +8.00	-14.00 to +8.00	-14.00 to +8.00	-14.00 to +8.00	-14.00 to +8.00	-14.00 to +8.00	
ADD POWER							
CYLINDER	to -7.00	to -7.00	to -7.00	to -7.00	to -7.00	to -7.00	
OTHER		Gray/Brown	Gray	Gray	Gray/Brown/G-15		
INVISIBLE MARKINGS	PC RLX	PT RLX	PX RLX	PV RLX	PPZ RLX	PD RLX	

OCCUPATION



PRODUCT FEATURES:

- ☐ FREEFORM PAL
- ☐ FREEFORM SV
- ☒ CONVENTIONAL PAL
- ☐ CONVENTIONAL SV
- ☒ EYEPOINT TECHNOLOGY®
- ☐ AS-WORN TECHNOLOGY®

NOTE:

- Recommended min fitting height is 16mm lower and 10mm upper (Base 2 not available in 2.25DP.)
- Cylinder to -6.00
- Prism to 4.00 Diopter
- Base Curve selection is based on the Total Near Rx vs. the Standard Distance Rx.

CR-39	TREATMENT				
	Clear	Transitions	XTRActive	Vantage	Polarized
TOTAL NEAR POWER RANGE	-6.00 to +6.00	--	--	--	--
DYNAMIC REDUCTION	-0.75, -1.25, -1.75, -2.25	--	--	--	--
BASE CURVE	2.5, 5, 7	--	--	--	--
OTHER/COATING	Hard Coated/Uncoated	--	--	--	--
INVISIBLE MARKINGS	-175	--	--	--	--

SuperLite™ 1.74	TREATMENT			
	Clear			
TOTAL NEAR POWER RANGE	--			
DYNAMIC REDUCTION	--			
BASE CURVE	--			
OTHER/COATING	--			
INVISIBLE MARKINGS	--			

SuperLite™ 1.60	TREATMENT				
	Clear	Transitions	XTRActive	Vantage	Polarized
TOTAL NEAR POWER RANGE	--	--	--	--	--
DYNAMIC REDUCTION	--	--	--	--	--
BASE CURVE	--	--	--	--	--
OTHER/COATING	--	--	--	--	--
INVISIBLE MARKINGS	--	--	--	--	--

DLC™ TRIVEX™	TREATMENT			
	Clear	Transitions	XTRActive	Polarized
TOTAL NEAR POWER RANGE	--	--	--	--
DYNAMIC REDUCTION	--	--	--	--
BASE CURVE	--	--	--	--
OTHER/COATING	--	--	--	--
INVISIBLE MARKINGS	--	--	--	--

SuperLite™ 1.67	TREATMENT				
	Clear	Transitions	XTRActive	Vantage	Polarized
TOTAL NEAR POWER RANGE	--	--	--	--	--
DYNAMIC REDUCTION	--	--	--	--	--
BASE CURVE	--	--	--	--	--
OTHER/COATING	--	--	--	--	--
INVISIBLE MARKINGS	--	--	--	--	--

POLYPLUS™ / POLYCARBONATE	TREATMENT				
	Clear	Transitions	XTRActive	Vantage	Polarized
TOTAL NEAR POWER RANGE	-3.00 to +6.00	--	--	--	--
DYNAMIC REDUCTION	-0.75, -1.25, -1.75, -2.25	--	--	--	--
BASE CURVE	5, 7	--	--	--	--
OTHER/COATING	Hard Coated: Non Primer:	--	--	--	--
INVISIBLE MARKINGS	P 175	--	--	--	--

Crizal AR Coating Hoya AR Coating Zeiss AR Coating iCoat AR Coating

All Shamir products are available with your lab's in-house AR Coatings. Please check with your lab to determine which specific manufacturer's AR coating is available with the Shamir lens and material that you are ordering, as a manufacturer may offer more than one AR coating.

EFFECTIVE DATE: SEPT. 1. 2012

OFF-GUI-REF-081612

